CE 2797

Conformité Européenne (European Conformity), MDR 2017/745
Conformité Européenne (Europäische Konformität), MDR 2017/745
Conformité Européenne (Conformité Européenne), MDR 2017/745
Conformité Européenne (Conformità europea), MDR 2017/745
Conformité Européenne (conformidad europea), MDR 2017/745
Conformité Européenne (europæisk overensstemmelse), MDR 2017/745
Conformité Européenne (eurooppalainen vaatimustenmukaisuus), MDR 2017/745
Conformité Européenne (Europos atitiktis), MDR 2017/745
Conformité Européenne (Europese conformiteit), MDR 2017/745
Conformité Européenne (Europeisk samsvar), MDR 2017/745
Conformité Européenne (Zgodność europejska), MDR 2017/745
Conformité Européenne (Conformidade Europeia), MDR 2017/745
Conformité Européenne (Conformitate europeană), MDR 2017/745
Conformité Européenne (i överensstämmelse med EU-direktiv), MDR 2017/745
Conformité Européenne (evropská konformita), MDR 2017/745
Conformité Européenne (európai megfelelés), MDR 2017/745

UK
CA 0086

United Kingdom Conformity Assessed

UKRP

United Kingdom Responsible Person:
Dot Medical Ltd. trading at Silk Point,
Queens Avenue, Macclesfield, SK10 2BB, United Kingdom

MAGiC
Magnetic Interventional Ablation Catheter


An Osypka product
Manufactured for:

Version: 2025-02-04

English	4
MAGiC Interventional Ablation Catheter Instructions for Use	
Deutsch	5
MAGiC Interventioneller Ablationskatheter Gebrauchsanweisung	
Français	7
MAGiC Cathéter d'ablation interventionnelle Instructions d'utilisation	
Italiano	8
MAGiC Catetere per ablazione interventistica Istruzioni per l'uso	
Español	10
MAGiC Catéter de ablación intervencionista Instrucciones de uso	
Dansk	11
MAGiC Interventionelt ablationskater Brugsanvisning	
Suomi	12
MAGiC Interventionaalinen ablaatiokateri Käyttöohjeet	
LIETUVIŲ K.	14
MAGiC Magnetinis intervencinės abliacijos kateteris Naudojimo instrukcijos	

Nederlands	15
MAGiC Interventionele ablatiekatheter Gebruiksaanwijzing	
Norsk	17
MAGiC Intervensjonelt ablasjonskater Bruksanvisning	
Polski	18
MAGiC Cewnik do ablacji interwencyjnej Instrukcja używania	
Português	20
MAGiC Cateter de ablação intervencionista Manual de instruções	
Română	21
MAGiC Cateter de ablație intervențională Instrucțiuni de utilizare	
Svenska	23
MAGiC Interventionell ablationskater Bruksanvisning	
Čeština	24
MAGiC Intervenční ablační katér Návod k použití	
Magyar	25
MAGiC Intervenciós ablációs katéter Használati utasítás	

	en • Date of manufacture de • Herstellungsdatum fr • Date de fabrication it • Data di costruzione es • Fecha de fabricación da • Fremstillingsdato fi • Valmistuspäivämäärä lt • Pagaminimo data nl • Fabricagedatum no • Dato for produksjon pl • Data produkcji pt • Data de fabrico ro • Data de fabricație sv • Tillverkningsdatum cs • Datum výroby hu • Gyártási dátum	LOT	en • Batch name de • Chargen-Bezeichnung fr • Désignation des charges it • Lotto, designazione es • Description of the load da • Lotbetegnelse fi • Valmistuserän kuvaus lt • Partijos pavadinimas nl • Lotaanduiding no • Partinavn pl • Nazwa partii pt • Designação do lote ro • Denumirea lotului sv • Partibeteckning cs • Číslo šarže hu • Gyártási tételszám	MD	en • Medical Device de • Medizinprodukt fr • Dispositif médical it • Dispositivo medico es • dispositivo médico da • Medicinsk udstyr fi • Lääkinnällinen laite lt • Medicinos prietaisas nl • Medisch hulpmiddel no • Medisinsk utstyr pl • Wyrób medyczny pt • Dispositivo Médico ro • Dispozitiv medical sv • Medicinteknisk produkt cs • Zdravotnický prostředek hu • Orvostechikai eszköz
	en • Manufacturer de • Hersteller fr • Fabricant it • Fabbicante es • Manufacturer da • Producent fi • Valmistaja lt • Gamintojas nl • Fabrikant no • Produsent pl • Producent pt • Fabricante ro • Producător sv • Tillverkare cs • Výrobce hu • Gyártó	REF	en • Order number de • Bestellnummer fr • Référence de commande it • Order number es • Reference number da • Bestillingsnummer fi • Tilausnumero lt • Užsakymo numeris nl • Bestelnummer no • Ordrenummer pl • Numer zamówienia pt • Número de encomenda ro • Numărul de comandă sv • Ordernummer cs • Číslo objednávky hu • Megrendelési szám	STERILEEO	en • Sterilisation with ethylene oxide de • Sterilisation mit Ethylenoxid fr • stérilisation à l'oxyde d'éthylène it • Sterilizzazione con ossido di etilene es • Esterilizado con óxido de etileno da • Steriliseret med ethylenoxidgas fi • Steriloitu etyleenioksidikaasulla lt • Sterilizavimas etileno oksidu nl • Gesteriliseerd met ethyleenoxidegas no • Sterilisering med etylenoksid pl • Wysterylizowane tlenkiem etylenu pt • Esterilizado com óxido de etileno ro • Sterilizare cu oxid de etilenă sv • Steriliserad med etylenoxidgas cs • Sterilizováno ethylenoxidovým plynem hu • Etilénoxid gázzal sterilizálva
	en • Consult Instructions for Use de • Gebrauchsanweisung beachten fr • Consulter les instructions d'utilisation it • Consult the instructions for use es • Consult the instructions for use da • Følg brugsanvisningen fi • Noudata käyttöohjetta lt • Žr. naudojimo instrukcijas nl • Gebruiksaanwijzing raadplegen no • Se bruksanvisningen pl • Przestrzegać instrukcji użytkowania pt • Observar as instruções de operação ro • Consultați instrucțiunile de utilizare sv • Beakta bruksanvisningen cs • Respektujte návod k použití hu • Tartsa be a használati utasítást		en • Contents de • Inhalt fr • Sommaire it • Contenuto es • Contenido da • Indhold fi • Sisältö lt • Turinys nl • Inhoud no • Innhold pl • Treść pt • Conteúdo ro • Conținut sv • Innehåll cs • Obsah hu • Tartalom		en • For one-time usage only de • Nur zum Einmalgebrauch fr • A usage unique it • Monouso es • De usar sólo una vez da • Kun til engangsbrug fi • Kertakäyttötuote lt • Tik vienkartiniam naudojimui nl • Uitsluitend bestemd voor eenmalig gebruik no • Kun til engangsbruk pl • Tylko do jednorazowego użytku pt • Utilizar apenas uma vez ro • Numai pentru o singură utilizare sv • Endast för engångsbruk cs • Pouze k jednorázovému použití hu • Csak egyszeri használatra
	en • Keep dry de • Trocken aufbewahren fr • Conserver au sec it • Mantenere secco es • Mantener seco da • Opbevares tørt fi • Suojaa märältä lt • Laikykite sausai nl • Tegen vocht beschermen no • Oppbevares tørt pl • Chronić przed wilgocią pt • Manter seco ro • Păstrați uscat sv • Förvaras torrt cs • Chraňte před mokrem hu • Szárazon tartandó!		en • Caution de • Vorsicht fr • Attention it • Attenzione es • Precaución da • Forsigtig fi • Varoitus lt • Įspėjimas nl • Let op no • Forsiktig pl • Ostrzeżenie pt • Cuidado ro • Atenție sv • Iakttag försiktighet cs • Obecný varovný hu • Figyelem!		en • Temperature limitation de • Temperaturbegrenzung fr • Limites de température it • Limiti di temperatura es • Limite de temperatura da • Temperaturbegrænsning fi • Lämpötilan rajoitus lt • Temperatūros apribojimas nl • Temperatuurgrens no • Temperaturbegrensning pl • Ograniczenie temperatury pt • Limite de temperatura ro • Limitarea temperaturii sv • Temperaturbegränsning cs • Teplotní omezení hu • Hőmérséklethatárok
	en • To be used before de • Verwendbar bis fr • Utiliser avant it • Utilizzabile fino a es • Utilizar antes da • Holdbar til fi • Viimeinen käyttöpäivä lt • Naudoiti prieš nl • Te gebruiken tot no • Brukes innen pl • Przydatne do użycia do pt • Utilizar até ro • A se utiliza înainte de sv • Används före cs • Spotřebuje do hu • Minőségét megőrzi		en • Do not use if package damaged de • Inhalt beschädigter Packung nicht verwenden fr • Ne pas utiliser si l'emballage est endommagé it • Non utilizzare se la confezione è danneggiata es • No usar si el paquete está dañado da • Indhold fra beskadiget emballage må ikke anvendes fi • Älä käytä vaurioituneen pakkauksen sisältöä lt • Nenaudokite, jei pakuotė pažeista nl • Inhoud van een beschadigde verpakking niet gebruiken no • Skal ikke brukes hvis emballasjen er skadet pl • Nie używać w przypadku uszkodzonego opakowania pt • Não usar o conteúdo de embalagem danificada ro • Nu utilizați dacă ambalajul este deteriorat sv • Använd inte innehållat vid skadad förpackning cs • Obsah poškozeného obalu nepoužívejte hu • Ha a csomagolás megrongálódott, a tartalmát ne használja fel		en • Do not sterilise anew de • Nicht erneut sterilisieren fr • Ne pas stériliser à nouveau it • Non sterilizzare di nuovo es • No esterilizar de nuevo da • Må ikke gensteriliseres fi • Ei saa steriloida uudestaan lt • Nesterilizuokite iš naujo nl • Niet opnieuw steriliseren no • Ikke steriliser på nytt pl • Nie sterylizować ponownie pt • Não esterilizar novamente ro • Nu sterilizați din nou sv • Får ej resteriliseras cs • Opakovaně nesterilizujte hu • Ne sterilizálja újból

	en • Single sterile barrier system with protective packaging inside de • Einfaches Sterilbarriersystem mit Schutzverpackung innen fr • Système de barrière stérile simple avec emballage de protection interne it • Sistema di barriera sterile singola con confezione protettiva interna es • Sistema de barrera estéril única con embalaje protector interior da • System med enkelt steril barriere med indvendig beskyttelseemballage fi • Yksinkertainen steriili estejärjestelmä ja sisäpuolinen suoja pakkaus lt • Vienkartinė sterili barjerinė sistema su apsaugine pakuote viduje nl • Toepassing van een enkelvoudige steriele barrière rondom een beschermende verpakking no • Enkelt sterilt barriersystem med beskyttende innvendig emballasje pl • System pojedynczej bariery sterylnej z wewnętrznym opakowaniem ochronnym pt • Sistema de barreira única estéril com embalagem de protecção interna ro • Sistem unic de barieră sterilă cu ambalaj de protecție în interior sv • System med enkel sterilbarriär och skyddsförpackning inuti cs • Jednoduchý systém sterilní bariéry s ochranným balením uvnitř hu • Egyszeres steril határolórendszer belső védőcsomagolással		en • Humidity limitation de • Feuchtigkeitsebegrenzung fr • Limitation de l'humidité it • Limitazione dell'umidità es • Limitación de la humedad da • Luftfugtighedsbegrænsning fi • Kosteusraja lt • Ribinė drėgmė nl • Vochtigheidslimiet no • Fuktighetsbegrensning pl • Zakres dopuszczalnej wilgotności pt • Limitação de humidade ro • Limitarea umidității sv • Fuktighetsbegränsning cs • Omezení vlhkosti hu • Páratartalom-korlátozás		en • Atmospheric pressure limitation de • Luftdruckbegrenzung fr • Limitation de la pression atmosphérique it • Limitazione della pressione atmosferica es • Limitación de la presión atmosférica da • Atmosfærisk trykbegrensning fi • Ilmanpaineraja lt • Ribinis atmosferos slėgis nl • Limiet van de luchtdruk no • Atmosfærisk trykkgrensning pl • Zakres dopuszczalnego ciśnienia atmosferycznego pt • Limitação de pressão atmosférica ro • Limitarea presiunii atmosferice sv • Atmosfärisk tryckbegränsning cs • Omezení atmosférického tlaku hu • Légköri nyomáskorlátozás
	en • Keep away from Sunlight de • Vor Sonnenlicht geschützt aufbewahren fr • Ne pas exposer aux rayons solaires it • Evitare l'esposizione alla luce del sole es • Evitar la exposición a los rayos solares da • Må ikke udsættes for sollys fi • Suojaa auringonvalolta lt • Saugokite nuo saulės spindulių nl • Tegen zonlicht beschermen no • Beskyttes mot sollys pl • Chronić przed nasłonecznieniem pt • Proteger da luz solar ro • Păstrați departe de lumina soarelui sv • Skyddas mot direkt solljus cs • Chraňte před slunečním světlem hu • Védjük a napfénytől!		en • Protect from Heat and Radioactive Sources de • Vor Hitze und radioaktiven Strahlenquellen schützen fr • Protéger de la chaleur et des sources radioactives it • Proteggere dal calore e dalle fonti radioattive es • Proteger del calor y las fuentes radiactivas da • Beskyt mod varme og radioaktive kilder fi • Suojaa lämmöltä ja radioaktiivisilta lähteiltä lt • Saugoti nuo karščio ir radioaktyvių šaltinių nl • Beschermen tegen hitte en radioactieve bronnen no • Beskytt mot varme og radioaktive kilder pl • Ochrona przed ciepłem i źródłami radioaktywnymi pt • Proteger do calor e de fontes radioativas ro • Protejați de căldură și surse radioactive sv • Skydd mot värme och radioaktiva källor cs • Chraňte před teplem a radioaktivními zdroji hu • Védelem a hőtől és a radioaktív forrásoktól		en • Medical prescription de • Verschreibungspflichtig fr • Prescription médicale it • Prescrizione medica es • Prescripción médica da • Kun på ordination fi • Lääkemääräys lt • Medicinos receptas nl • Medisch voorschrift no • Medisinsk resept pl • Recepta lekarska pt • Prescrição médica ro • Prescripție medicală sv • Recept cs • Lékařský předpis hu • Orvosi vény
	en • Unique Device Identifier de • Eindeutige Produktidentifizierung fr • Identificateur unique du dispositif it • Identificatore univoco del dispositivo es • Identificador único del dispositivo da • Unik identifikator for enheden fi • Yksilöllinen laitetunniste lt • Unikalus įrenginio identifikatorius nl • Unieke apparaatidentificatiecode no • Unique Device Identifier pl • Unikalny identyfikator urządzenia pt • Identificador único do dispositivo ro • Identificatorul unic al dispozitivului sv • Unik enhetskod för identifiering av enheten cs • Jedinečný identifikátor zařízení hu • Egyedi eszközazonosító		en • United Kingdom Conformity Assessment de • Vereinigtes Königreich Konformitätsbewertung fr • Royaume-Uni Évaluation de la conformité it • Regno Unito Valutazione della conformità es • Reino Unido Evaluación de la conformidad da • Storbritannien Overensstemmelsesvurdering fi • Yhdistynyt kuningaskunta Vaatimustenmukaisuuden arviointi lt • Jungtinės Karalystės atitikties nl • Verenigd Koninkrijk Conformiteitsbeoordeling no • Samsvarsvurdering for Storbritannia pl • Wielka Brytania Ocena zgodności pt • Avaliação da conformidade no Reino Unido ro • Regatul Unit Evaluarea conformității sv • Förenade kungariket Bedömning av överensstämmelse cs • Spojené království Posuzování shody hu • Egyesült Királyság Megfelelőségértékelés		

Operating instructions

MAGiC™

Magnetic Interventional Ablation Catheter

English

Caution: Please read these instructions carefully before using the product.

UDI: Basic UDI-DI: 4044508MAGICPM
SSCP: The Summary of Safety and Performance is publicly available at <https://ec.europa.eu/tools/eudamed>
If not accessible in EUDAMED, the information shall be made available to the public upon request
eIFU: <http://www.stereotaxis.com/referencelibrary>

Components
Magnetic Interventional Ablation Catheter (MAGiC™)

1. Device Description
MAGiC is a sterile, magnetically guided mapping and ablation catheter. MAGiC is intended to provide sensing and mapping of the electrical signals in the heart for diagnostic purposes. It has four electrodes for pacing and sensing with an open irrigation tip electrode for delivering ablative energy to the cardiac tissue. MAGiC is intended to be used with a compatible Stereotaxis Robotic Magnetic Navigation System (RMNS) and the Cardiodrive® catheter advancement system (CAS).

2. Essential use information

Characteristic	Specification	Accuracy
Catheter shaft diameter	8.5 F	2.67-2.87 mm
Catheter usable length	135 cm	133 cm ± 2 cm
Tip electrode length	3.5 mm	3.5 mm ± 0.1 mm
Tip electrode material	Gold	
Electrode material (except tip)	Platinum/Iridium	
Ring electrode length	1.5 mm	1.5 mm ± 0.1 mm
Electrode Spacing	2-3-1 mm	Tolerance: +0.3/-0.2 mm

3. Combinations with other medical devices

The following devices are compatible with MAGiC:

Legal Manufacturer	Device
Livetec Ingenieurbuero GmbH	HAT 500® RF Generator, HAT 500 Irrigation Pump, Remote Control, tubing
Abbott	Ampere™ RF Generator, compatible irrigation pump, tubing
Boston Scientific	Maestro Cardiac Ablation System, compatible irrigation pump, tubing
Biosense Webster	SMARTABLATE™ System (RF Generator, pump, and remote control), tubing
Stereotaxis	Robotic Magnetic Navigation System (RMNS)
Stereotaxis	Cardiodrive® and QuikCAS™
OSYPKA AG	iCONNECT™ and iCONNECT Cables

4. Indications for use

MAGiC is intended for use in cardiac electrophysiology for mapping, delivering diagnostic pacing stimuli, and for patients indicated for RF catheter ablation of cardiac arrhythmias. The indication for treatment of these arrhythmias is described in the guidelines and consensus statements from major cardiology societies, including the European Society of Cardiology.

5. Contraindications

The catheter should not be used:

- in patients with intracardiac mural thrombus;
- in patients where the catheter would need to cross a prosthetic valve;
- in the coronary arteries;
- with patients having a systemic infection;
- in patients with history of sensitivity to foreign objects or extreme allergies;
- in patients with histological or anatomical abnormalities that may lead to post-operative complications such as bleeding diathesis, diminished resistance to infection;
- in patients with hemodynamic instability;
- with patients unable to receive heparin or an acceptable alternative to achieve adequate anticoagulation.

6. Intended purpose

MAGiC together with Stereotaxis RMNS is intended to be used during radio-frequency ablation to produce local lesions in cardiac tissue for the treatment of cardiac arrhythmias, for electrophysiological mapping, and for delivering stimulation pulses.

7. Expected clinical benefit

Ablation of the tissues that cause or support an arrythmia has been shown in many patients to eliminate the arrythmia and return the heart to a normal sinus rhythm.

8. Patient population

MAGiC is intended to be used with patients that have cardiac arrythmia whose health status is suitable for a cardiac ablation procedure and do not have any of the listed contraindications.

9. How supplied

The products are sterilized with ethylene oxide. Packaging is designed to maintain sterility according to the expiration date on the label. Do not use if the expiration date has passed, the package is opened or damaged, or if the labelling is incomplete or illegible.

10. Handling and storage

Store at ambient temperatures in a dry and dark location (10-25°C, 30-60%RH, 60-120kPa). Keep dry and protect from direct sunlight. Do not expose to generated sources of gamma, neutron, α -particle and β -particle radiation.

11. Reuse precaution statement

Not autoclavable. Single use only. Do not re-sterilize product. Re-use bears the risk of infection and the possibility of malfunction.

12. Warnings

- **For single use only.** Do not reuse, reprocess, or resterilize. Reuse can compromise the catheter performance characteristics and may result in patient injury or infection.
- This catheter should be used only by physicians trained in electrophysiology techniques.
- This catheter should be used only in procedures that use the Stereotaxis RMNS.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.
- Use of devices or accessories other than those specified could result in improper operation.
- Proper sterile technique must be used when removing the catheter from the packaging and during operation of the catheter to prevent exposing the patient to the possibility of infection.
- Do not modify the product in any way.
- Do not advance through a prosthetic valve.
- Do not use after the expiry date on the package label.

13. Precautions

- The catheter is supplied STERILE in an unopened package. Verify package integrity has been maintained and the sterility of the device has not been compromised.
- Under sterile conditions, examine the device for defects and verify proper device function and integrity.
- If damage to the sterile barrier or device are found, do not use the device. Contact your Stereotaxis representative and exchange the device for a new one.
- If higher power and times are used than indicated in the RF parameters table, there is a possibility of steam pop occurring.
- Always use an irrigation pump with the catheter.
- Always use the RMNS, Cardiodrive CAS, and Cardiodrive Sterile Components with the catheter.
- Monitor the amount of irrigation fluid used to prevent fluid overload in the patient.
- Do not expose the catheter to organic solvents.
- Do not autoclave the catheter.
- Do not immerse the electrical connections in water.
- Do not clean the handle with liquids.
- Ensure that all air is removed from pump tubing and the catheter through proper flushing.
- Use intravenous heparin to avoid thromboemboli, particularly during left-sided heart access.
- Always maintain a base irrigation flow of 2ml/min.
- Do not use in or near MRI equipment.
- Inspect irrigation saline for bubbles. Air bubbles may cause emboli.
- Reduce the magnetic field when inserting or removing the catheter to prevent damage to the catheter.
- Ablating near the SA node or AV node may cause damage requiring temporary or permanent pacing.
- Use fluoroscopic imaging and electrogram signals to monitor catheter movements in order to reduce the risk of tissue injury.
- Before use, check that the irrigation ports are open and irrigation fluid is flowing out of the holes by flushing with saline.
- Do not exceed 24 hour use of the catheter.
- Observe error messages on the ablation generator and consult the IFU for guidance.
- Follow the respective operating instructions when connecting and using external devices.

14. Intended Users (=Use and handling)

The product may only be used in a medical treatment facility that has been specifically set up for the appropriate application and by trained staff (electrophysiology cath lab specialists). Physicians and staff using the MAGiC catheter must be trained in the use of the Stereotaxis Magnetic Navigation System. It is the responsibility of the physician to select the medically appropriate procedure for the patient. These operating instructions are intended to provide general information on handling the product. The product must not be modified in any way. Warnings, general information, and precautionary measures must be observed. Failure to observe these instructions can result in procedural complications or malfunctioning of or damage to the device.

15. Preparation for use

1. Carefully remove the catheter from the package using sterile techniques and place in a sterile work area.
2. Inspect the catheter prior to use to confirm it is undamaged. If damage is detected (e.g., bends, kinks, etc.), do not attempt to repair the catheter. Use a new catheter for the procedure and return the damaged catheter to Stereotaxis, Inc.
3. Remove a compatible introducer using sterile techniques and place in a sterile work area.

16. Instructions for use

1. Using aseptic techniques create vascular access in a large central vessel.
2. Prepare and flush the Irrigation Tubing Set according to the device IFU.
3. Connect the Irrigation Tubing Set to the catheter using standard luer fittings.
4. Flush the catheter and tubing to ensure the removal of any trapped air bubbles and to verify that the irrigation holes are open and irrigation fluid is flowing out of the holes at a flow rate of 2 ml/min.
5. Confirm the magnetic fields are removed.
6. Insert the catheter through the QuikCAS hemostasis adapter, following the QuikCAS Instructions for Use.
7. Insert the tip of the catheter into the hemostasis valve of the Introducer.
8. Lock the hub of the QuikCAS hemostasis adapter to the hub of the introducer.
9. Slide the insertion tube of the QuikCAS hemostasis adapter into the hemostasis valve of the introducer to open the hemostasis valve.
10. Advance the catheter until the soft distal section with the magnets has been inserted into the introducer (about 10cm).
11. Aspirate the introducer to remove any air that may have entered the introducer with the catheter's insertion.
12. Retract the handles of the QuikCAS hemostasis adapter to withdraw the insertion tube from the hemostasis valve.
13. Advance the catheter into the heart using fluoroscopic guidance.
14. Connect the rear disk of the QuikCAS hemostasis adapter to the QuikCAS.
15. Place the catheter into the QuikCAS and lock the catheter into the drive.
16. Using robotic magnetic navigation, perform the necessary mapping, sensing, and pacing required to diagnose the arrhythmia and identify the location for therapy delivery.

17. Deliver RF energy as necessary to treat the arrhythmia.
18. Reduce the magnetic field using the Apply magnetic field/Reduce magnetic field icon on the Main toolbar of the Navigant™ Workstation Software or retract the RMNS prior to removing the catheter.
19. When removing the catheter (during or after the procedure), remove the hemostasis adapter from the QuikCAS. Gently withdraw the catheter following best practices.

17. RF Ablation

Connect the catheter to the appropriate input on the iCONNECT interface box or the RF generator and use a compatible cable. An indifferent electrode patch must be placed on the patient and connected to the RF generator. Upon inserting the catheter, ensure that the displayed impedance is in an acceptable range and that the displayed temperature is near body temperature (~37°C). If impedance is not acceptable, the generator will not operate. Consult the applicable IFU for questions about the RF generator.

18. RF Parameters

Characteristic	Specification
Power Range	15-50W
Temperature Monitoring	<50°C
Application time	5-60 seconds for power of 30W or less 5-20 seconds for power greater than 30W
Irrigation Rate	10 ml/min

19. Potential adverse events

Adverse events may result from proper and improper use of the device and may also be due to the procedure itself. Follow instructions for use carefully. Potential adverse events related to the device and/or procedure vary in frequency and severity up to and including Death and include, but are not limited to:

atrioesophageal fistula	embolism	infection	pulmonary vein stenosis
bleeding	femoral arteriovenous fistula	inguinal hematoma	stroke
cardiac arrhythmia	femoral bleeding	myocardial infarction	thrombosis
cardiac tamponade	heart failure	perforation	transient phrenic palsy
conduction block	hemothorax	pericardial effusion	
diaphragm palsy	inflammation	phrenic nerve injury	

20. Disposal

The product may be contaminated after use and should be considered a potential biohazard. Handle and dispose-of in accordance with medical practice and applicable local, state, and federal laws and regulations.

21. Trademarks

Stereotaxis, the Stereotaxis logo, MAGiC, Cardiodrive, QuikCAS, Genesis Robotic Magnetic Navigation System, and Niobe Magnetic Navigation System are trademarks or registered trademarks of Stereotaxis, Inc. in the USA and other countries. All brand names, product names, or trademarks are the property of their respective owners.

22. Stereotaxis technical support

For technical support, please contact Stereotaxis TeleRobotic Support Team (TST) at 1-866-269-5268 or 1-314-678-6200 or email tst@stereotaxis.com.

23. Notice to the user and/or patient

Any adverse event or malfunction that occurred in relation to the device should be reported to the manufacturer. Serious adverse events should be reported to the competent authority of the Member State in which the user and/or patient is established.

24. Disclaimer of warranty and limitation of liability

STEREOTAXIS MAKES NO WARRANTIES, EXPRESS OR IMPLIED, WITH REGARD TO THE DEVICE DESCRIBED IN THIS DOCUMENT. STEREOTAXIS DISCLAIMS ALL REPRESENTATIONS OR WARRANTIES, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR USE, TITLE OR NON-INFRINGEMENT, ARISING BY STATUTE OR IN LAW, OR ARISING FROM A COURSE OF CONDUCT, COURSE OF DEALING OR USAGE OF TRADE. STEREOTAXIS, INC. SHALL NOT BE LIABLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. STEREOTAXIS IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES RESULTING FROM REUSE OF THIS DEVICE.

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