MAGiC SWEEP™

INSTRUCTIONS FOR USE

REF 200871-001


English	4	Lietuvių k.	13
中文	5	Magyar	14
繁體中文	5	Nederlands	15
Česky	6	Polski	16
Dansk	7	Português	17
Deutsch	9	Русский	19
Español	10	Română	20
Français	11	Suomi	21
Italiano	12	Svenska	22

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

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
DEVICE INTRODUCTION

The MAGiC Sweep™ EP Mapping Catheter is a flexible, insulated radiopaque catheter constructed of thermoplastic elastomer material with a variable number of noble metal electrodes. The tip of the steerable catheter is magnetically enabled and manipulated by the Stereotaxis Robotic Magnetic Navigation System. It can be used for recording intracardiac signals and for cardiac stimulation during electrophysiology studies or to map regions of heart.

INDICATIONS

The MAGiC Sweep™ EP Mapping Catheter is intended for intracardiac electrophysiological recording and/or stimulation for pacing in the heart.

PATIENT TARGET GROUP

The MAGiC Sweep™ EP Mapping Catheter targets patients with cardiac arrhythmia in need of cardiac stimulation or mapping regions of their hearts.

INTENDED USERS

The MAGiC Sweep™ EP Mapping Catheters intended users are physicians trained in the technique of angiography, electrophysiology, intracardiac recording and stimulation, and in the use of the Stereotaxis Robotic Magnetic Navigation System.

CONTRAINDICATIONS

The MAGiC Sweep™ EP Mapping Catheters are contraindicated for:

- Electrophysiology studies are contraindicated when acute factors make the findings unrepresentative of the patient's usual state (e.g. electrolyte abnormality, acute ischemia, drug toxicity)
- When a patient's underlying cardiac disease makes it likely that induced arrhythmias will be extremely difficult to terminate and carry a high risk of death (e.g. acute myocardial infarction, unstable angina, hemodynamic instability)

WARNINGS and PRECAUTIONS

- This device should only be used by physicians trained in the techniques of angiography, electrophysiology, and intracardiac recording and stimulation
- Read the directions prior to use
- Do not alter this device
- United States Federal law restricts this device to sale by or on the order of a physician
- This device should only be used with equipment that complies with international safety standards
- Do not use this catheter for intracardiac ablation or for DC ablation.
- Do not use for internal cardioversion
- Do not force the catheter through the vessel
- Catheter advancement should be done under fluoroscopic guidance.
- This device is intended for one time use only. Do not reprocess or reuse. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another
- Store in a cool, dry place
- Follow all warnings and precautions applicable to the Stereotaxis Robotic Magnetic Navigation System, especially those regarding magnetic objects in the catheter laboratory. Refer to the Robotic Magnetic Navigation System Users Manual for user precautions
- Always reduce the magnetic fields of the Robotic Magnetic Navigation System before attempting to remove the magnetic catheter
- This device contains small permanent magnets. Keep six inches away from magnetically susceptible devices in order to avoid any potential interaction

CLINICAL BENEFIT

The MAGiC Sweep™ EP Mapping Catheter clinical benefits include capturing critical cardiac information to deliver a more precise therapy or treatment to fix arrhythmias.

TRANSPORTATION AND STORAGE

- Temperature: -29°C ~ 55°C
- Relative Humidity: Below 85%
- Keep away from light

CATHETER HANDLING AND CARE

- Personnel handling the electrophysiology catheter should follow appropriate sterile techniques
- Excessive bending or kinking of the catheter may cause damage to the catheter
- Do not allow moisture onto the connectors on the catheter
- Do not wipe the catheter with organic solvents, such as alcohol

- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy

USE BY DATE

Use the device prior to the "Use By" date on the package label.

DISPOSAL

Recycle components or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

ADVERSE EVENTS

Potential adverse events associated with cardiac catheter procedures include, but are not limited to the following:

- The risks of use of electrophysiology catheters include those risks related to cardiac catheterization in general, such as thromboembolisms, cardiac perforation, tamponade and infection. The induction of atrial fibrillation, ventricular tachycardia (VT) requiring cardioversion, and ventricular fibrillation (VF) can be risks associated with electrical stimulation
- Radiation exposure: electrophysiologic studies procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging
- Vascular access: (e.g. bleeding, hematomas, deep venous thrombosis, arterial or cardiac tissue perforation, cardiac tamponade, pulmonary embolism, thromboembolism, myocardial infarction, stroke, arteriovenous fistula, pneumothorax, local or systemic infection, and death)
- Catheter manipulation: (e.g. valvular damage, micro emboli, perforation of the coronary sinus or myocardial wall, coronary artery dissection, thrombosis)

Any serious incident that has occurred in relation to the device should be reported to Access Point Technologies EP and the competent authority of the Member State the patient/user is established.

INSTRUCTIONS FOR USE

1. Inspect the catheter packaging before opening. Do not use if the package is open or damaged.
2. Use aseptic technique to remove the catheter from the package and place it on a sterile working area.
3. Carefully inspect the electrode integrity and overall condition of the catheter. Do not use the catheter if the electrodes or tip are loose, irregular or otherwise visibly damaged.
4. Perform vascular access utilizing a sterile technique. The catheter may be used in the femoral, brachial, subclavian, or jugular access sites.
5. Connect the catheter to a compatible extension cable by aligning the marks of the catheter's connector with the corresponding marks on the cable. Press the connectors together. Do not force the connection.
6. Connect the lead pins of the Catheter Connection Cable to the appropriate electronic recording and/or stimulation equipment.
7. Insert the catheter through the QuikCAS™ hemostasis adapter, following the QuikCAS Instructions for Use.
8. Insert the tip of the catheter into the hemostasis valve of the Introducer.
9. Lock the hub of the QuikCAS hemostasis adapter to the hub of the introducer and slide the insertion into the hemostasis valve of the introducer to open the hemostasis valve.
10. Advance the catheter until the soft distal section with the magnets has been inserted into the introducer (about 10 cm).
11. Aspirate the introducer to remove any air that may have entered the introducer with the catheter's insertion.
12. Retract the handles of the QuikCAS hemostasis adapter to withdraw the insertion tube from the hemostasis valve.
13. Advance the catheter into the heart using fluoroscopic guidance.
14. Connect the rear disk of the hemostasis adapter to the QuikCAS.
15. Place the catheter into the QuikCAS and lock the catheter into the drive.
16. Advance the catheter into the desired area of the heart using fluoroscopic and ECG guidance. Using robotic magnetic navigation, perform the mapping and sensing required to diagnose the arrhythmia and identify locations for therapy delivery
17. Reduce the magnetic field or retract the Robotic Magnetic Navigation System prior to removing the catheter.
18. When removing the catheter (during or after the procedure), remove the hemostasis adapter from the QuikCAS. Gently withdraw the catheter.
19. Do not re-sterilize or reuse.
20. Dispose of the catheter according to hospital biohazard requirements.

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Access Point Technologies EP Inc.