Map-iT[®] Mapping and Ablation Catheter Extension Cable

INSTRUCTIONS FOR USE

Issue Date: 1/2026

Device Introduction

The Map-iT[®] Mapping and Ablation Catheter Extension Cable is an extension cable for the mapping and ablation catheter. It is comprised of an elongated and tubular cable body, two adapter connectors with different numbers of pins.

The Map-iT[®] Mapping and Ablation Catheter Extension Cable is designed for connecting the catheter to the RF generator or the electrocardiogram measuring instrument, and for electrical signal transmission during electro-physiology studies including mapping, pacing, and radiofrequency ablation.

The Map-iT[®] Ablation Catheter Extension cable is intended for use with a Map-iT[®] Ablation Catheter and compatible RF generator with the following specifications.

-Frequency: 420-575Kilohertz

-Voltage: 0-220Volts

-Current: 0-2.0Amps

-Power: 150Watts Maximum

-Peak value: 100Volts (40-60Volts is normal)

-Accessory current: 0.2-0.6Amps

These output ranges represent the device capability, not the actual settings. The actual output settings will be determined by physician experience and based on clinical requirements. APT EP recommends operating the compatible RF generator in the Temperature control mode when using systems which are capable of delivering high power.

Normal conditions of use:

Operating temperature: 10° C ~40° C (50° ~140° F), Humidity rate: ≤ 85%, Pressure: Normal atmospheric pressure.

Indication

The Map-iT[®] Mapping and Ablation Catheter Extension Cable is indicated for connecting the mapping or ablation catheter to the multichannel physiological recorder or the RF generator. It is used for transporting electrical signals during diagnostic electrophysiology studies or cardiac ablation procedures.

Contraindications

There are no known contraindications for this cable.

Warnings

-  Store in a cool, dark, dry place
-  Do not alter this device
-  Inspect components prior to use
-  Do not reuse this device. Thorough cleaning of biological and foreign material is not possible. Mechanical damage and /or adverse patient reactions may result from the reuse of this device.
-  Use of this device should only extend to those physicians who are skilled in the techniques of transvenous intracardiac studies, temporary pacing and RF ablation. The RF generator, accessories and electrode catheters should be tested before clinical operator.
-  Proper functioning of the Map-iT[®] Mapping and Ablation Catheter Cable requires that you prevent moisture from contacting the connector
-  Only Map-iT[®] Mapping and Ablation Catheters can be used with the Mapping and Ablation Catheter Extension Cable.
-  Map-iT[®] Ablation Catheter Extension Cables are intended for use with a compatible RF generator. The RF generator should comply with 59.103.5 of IEC 60601-2-2:2006/ED4. Refer to the RF generator manual for detailed generator operation instructions for RF ablation (the RF generator catalog list: Medtronic Atakr RF generator, Johnson & Johnson Stockert RF generator, IBI RF generator, EPT RF generator, etc).
-  These instructions are general in nature. The physician may wish to defer from these instructions in light of personal clinical experience.
-  After use, dispose of the product and packaging in accordance with hospital, administrative and/or local government policies.

SUGGESTED DIRECTIONS FOR USE

1. Use this device to combine a compatible RF generator and appropriate ablation catheter, or to combine a multichannel physiological recorder and appropriate mapping catheter.
2. The connecting devices for the Map-iT[®] Mapping and Ablation Catheter Extension Cables include: radiofrequency (RF) generator, multi-channel physiological recorder, pacing equipment and cardiac programmable stimulator.
3. To record intracardiac electrograms, connect the extension cable (s) to the catheter. Observe the polarity of the proximal end connector pins of the extension cable (s) when connecting to an ECG amplifier that complies with international safety standards.
4. Consult the compatible RF generator instructions for use and multi-channel physiological recorder instructions for use for the proper connection of the patient grounding patch
5. Use care to isolate any unused connector pins. This will reduce the chances of developing accidental current pathways to the heart.

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Access Point Technologies EP Inc.