Map-iT[®]

Diagnostic Mapping Catheter

INSTRUCTIONS FOR USE


Access Point Technologies EP, Inc.


Issue Date: 09/2025

Device Introduction:

The Map-iT[®] Diagnostic Mapping Catheters are flexible, insulated radiopaque catheters constructed of thermoplastic elastomer material with a variable number of noble metal electrodes. The tip of the steerable catheters may be manipulated with the control mechanism located in the handle at the proximal end of the catheter. It can be used for recording intracardiac signals and for cardiac stimulation during electrophysiology studies or to map regions of heart.

Fixed Catheter

Steerable Catheter

Indication:

The Map-iT[®] Diagnostic Mapping Catheters can be used in the evaluation of a variety of cardiac arrhythmias from endocardial and intravascular sites.

Contraindications:

The Map-iT[®] Diagnostic Mapping Catheters are contraindicated for:

- Electrophysiology studies are contraindicated when acute factors make the findings unrepresentative of the patient's usual state (e.g. electrolyte abnormality, acute ischemia, drug toxicity).
- When the patient's underlying cardiac disease makes it likely that induced arrhythmias will be extremely difficult to terminate and carry a high risk of death (e.g. acute myocardial infarction, unstable angina, hemodynamic instability).
- Pregnant patients
- Patients with active systemic infection
- Patients with ventriculotomy or atriotomy recently
- Patients who are uncooperative

Warnings and Precautions:

- This device should only be used by physicians thoroughly trained in the techniques of angiography, electrophysiology, and intracardiac recording and stimulation
- Read the directions prior to use.
- Do not alter this device.
- United States Federal law restricts this device to sale by or on the order of a physician.
- This device should only be used with equipment that complies with international safety standards.
- Do not use this catheter for intracardiac ablation or for DC ablation.
- Do not use for internal cardioversion.
- If defibrillation is needed, disconnect the catheter.
- Do not force the catheter through the vessel.
- Catheter advancement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- This device is intended for one time use only. Do not reprocess or reuse. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another.
- Store in a cool, dry place.
 - Temperature: -29°C ~ 55°C
 - Relative Humidity: Below 85%
 - Keep away from light
- Catheter handling and care:
 - DO NOT insert or withdraw the catheter without straightening the catheter tip, as confirmed via fluoroscopy.
 - Do not deflect the catheter tip while it is constrained within a sheath.
 - Personnel handling the electrophysiology catheter should wear gloves.
 - Excessive bending or kinking of the catheter may cause damage to the catheter. Manual pre-bending of the distal curve can damage the steering mechanism and may cause patient injury.
 - Do not allow moisture onto the connectors on the catheter.
 - Do not wipe the catheter with organic solvents, such as alcohol.

- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Use by Date

Use the device prior to the “Use By” date on the package label.

Disposal

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

Adverse Events

Potential adverse events associated with cardiac catheter procedures include, but are not limited to the following:

- The risks of use of electrophysiology catheters include those risks related to cardiac catheterization in general, such as thromboembolisms, cardiac perforation, tamponade and infection. The induction of atrial fibrillation, ventricular tachycardia (VT) requiring cardioversion, and ventricular fibrillation (VF) can be risks associated with electrical stimulation.
- Radiation exposure: electrophysiologic studies procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging.
- Vascular access: (e.g. bleeding, hematomas, deep venous thrombosis, arterial perforation, cardiac tissue perforation, cardiac tamponade, pulmonary embolism, thromboembolism, myocardial infarction, stroke, arteriovenous fistula, pneumothorax, local or systemic infection, and death).
- Catheter manipulation: (e.g. valvular damage, micro emboli, perforation of the coronary sinus or myocardial wall, coronary artery dissection, thrombosis).

Instructions for Use

General

1. Inspect the catheter packaging before opening. Do not use if the package is open or damaged.
2. Use aseptic technique to remove the catheter from the package and place it on a sterile working area.
3. Carefully inspect the electrode integrity and overall condition of the catheter. Do not use the catheter if the electrodes or tip are loose, irregular or otherwise visibly damaged.
4. Perform vascular access utilizing a sterile technique. The catheter may be used in the femoral, brachial, subclavian, or jugular access sites.
5. Connect the catheter to a compatible extension cable.
 - a. Connect the cable by aligning the marks of the catheter's connector with the corresponding marks on the cable.
 - b. Press the connectors together. Do not force the connection.
6. Connect the lead pins of the Catheter Connection Cable to the appropriate electronic recording and/or stimulation equipment.
7. Advance the catheter into the desired area of the heart using fluoroscopic and ECG guidance.

Note: Steerable Handle Control Operation

 - a. Push the Tip Deflection Plunger Knob forward to deflect the tip of the catheter. The tip can be deflected to achieve the labeled curve shape, pull back on the Tip Deflection Plunger Knob.
 - b. Before withdrawing the catheter, fully pull back on the Tip Deflection Plunger Knob to un-deflect the tip of the catheter and verify by fluoroscopy that the tip is in a straight position.
8. Do not re-sterilize or reuse.
9. Dispose of the catheter according to hospital biohazard requirements.

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Access Point Technologies EP Inc.