Map-iT[®]

Mapping Catheter Extension Cable

INSTRUCTIONS FOR USE

Issue Date: 09/2025

INTRODUCTION

The Map-iT® Diagnostic Mapping Catheter Extension Cable is an extension designed for connecting the catheter to an electrocardiogram recording system and for electrical signal transmission during electrophysiology studies including mapping and pacing. It is comprised of an elongated and tubular cable body, electrical connectors on each end with the number of pins as specified by the model.

INDICATIONS FOR USE

The Map-iT® Diagnostic Mapping Catheter Extension Cable is indicated for connecting the mapping catheter to the multichannel physiological recorder. It is used for transporting electrical signals during diagnostic electrophysiology studies. This cable may be reprocessed up to 10 times subject to the cleaning and sterilization restrictions herein.

CONTRAINDICATIONS

There are no known contraindications for this cable.

Warnings and Precautions

- ⚠ Use the device prior to the “Use By” date on the package label.
- ⚠ Inspect cable for damage prior to use. Do not use if the cable is damaged (i.e. cuts, kinks, nicks, or crushed or elongated sections).
- ⚠ Do not use if a break occurs in the cable wire or the cable becomes otherwise discontinuous, arcing may occur and may burn the patient or create a fire.
- ⚠ If used in the presence of electrical equipment, noise may be induced into signals conducted by the cable.
- ⚠ Only Map-iT® Diagnostic Mapping Catheters can be used with the Map-iT® Diagnostic Mapping Catheter Extension Cable.
- ⚠ Do not alter this device.
- ⚠ Clean and re-sterilize the cables per the instructions provided.
- ⚠ Failure to properly clean cables could lead to inadequate re-sterilization.
- ⚠ Automated cleaning of the cables is not recommended and may damage the cable.
- ⚠ Do not immerse the cable connectors in fluid.
- ⚠ Do not expose the cable to organic solvents.
- ⚠ Do not reuse the cable after eleven (11) uses. Discard the cable after eleven (11) uses. Exceeding ten (10) re-sterilization cycles may impair the proper functioning of the cable.
- ⚠ Patient or operator injury can result from improper handling or storage.

PRODUCT CONDITION AND STERILIZATION

This cable is provided sterile (by EO) and may be re-sterilized for up to 10 times, per the instructions below, using H2O2 sterilization.

CLEANING AND RE-STERILIZATION

Cables should be cleaned and re-sterilized as soon as is practical after use. The cable should be cleaned and re-sterilized using the method specified in Table 1.

After cleaning, cables should be packaged in a breathable pouch prior to re-sterilization. Cables should be coiled into a loop when placing into the pouch. Use a pouch large enough to prevent stress to the pouch seams or excessive bending of the cable.

Sterilization cycles should be performed using techniques and parameters compliant with the sterilization equipment manufacturer's instructions and facility protocols.

The cleaning and re-sterilization parameters provided in Table 1 have been validated by the manufacturer.

Table 1

CLEANING AND STERILIZATION INSTRUCTIONS	
1.	Point-of-Use (in suite) a. Do not disassemble device. b. Remove gross soil and residual soil from surfaces and edges using a low-lint wipe moistened with water.
2.	Containment and Transportation a. Devices should be cleaned within 30 minutes of use to prevent drying prior to cleaning. Do not allow contaminated devices to dry prior to reprocessing. b. Universal precautions should be observed by all hospital personnel exposed to contaminated or potentially contaminated devices.
3.	Perform Device Initial Wipe a. Wipe cable with a moistened soft-bristled nylon brush (such as an appropriately sized cannula brush) moistened with utility water, low-lint wipes, and/or gloved hands to aid in the removal of visible gross soil. b. Inspect device for damage or malfunction. c. Be careful to avoid getting liquid inside the electrical connectors.

4.	Perform Device Scrub:
a.	Prepare a solution with neutral detergent and warm tap water following the detergent manufacturer's recommendation. The recommended detergent is Polystyca 2x Neutral mixed with warm tap water (<55°C) at 1 mL/L concentration. Note: Alkaline or acid based detergents may harm the device and should not be used.
b.	Utilizing the prepared neutral detergent solution, scrub entire surface area of device for a minimum of 2 minutes using a soft nylon instrument brush and low-lint wipes; give special attention to challenging areas using an appropriately sized channel brush.
c.	Be careful to avoid getting liquid inside the electrical connectors.
d.	Following the scrub, inspect the device for visible soil residue.
e.	If visible soil residue is observed on the device, repeat operation 3 step a. and operation 4 steps a. through c. until there is no visible soil residue.
5.	Perform Device Final Wipe:
a.	Wipe device for a minimum of 2 minutes with a low-lint wipe, moistened with water, paying special attention to challenging areas.
b.	Inspect for visible soil residue.
c.	If visible soil residue is observed on the device, repeat operation 3 and Operation 4 steps a. through c. until there is no visible residue.
6.	Dry Device
a.	Dry device using compressed air and a dry lint-free cloth until the device exterior surfaces are visually dry; articulate cable while drying to ensure all moisture is removed.
7.	Package and label device for sterilization
a.	Coil each cable so that diameter of coil is between 8-11 inches.
b.	Use a Tyvek/Polyethylene pouch (i.e. Vis-U-All Low Temperature Sterilization Pouches) large enough to prevent stress to the pouch seams or excessive bending of the cable.
c.	Place cable into pouch so that no excessive bending is required.
d.	Seal pouch
e.	Label pouch per institutional protocols/procedures.
8.	Sterilize
a.	Amsco V-PRO maX Low Temperature Sterilization System
i.	Load sterilizer per sterilization equipment operator instructions.
ii.	Select pre-programmed cycle: Non-Lumen Cycle (approximately 28 minutes to complete) using VAPROX HC Sterilant in specially designed disposable cartridges.
iii.	Operate machine per sterilization equipment operator instructions.
b.	STERRAD NX Sterilization System
i.	Load sterilizer per sterilization equipment operator instructions.
ii.	Select cycle: Standard
iii.	Operate machine per sterilization equipment operator instructions.
c.	STERRAD 100NX Sterilization System
i.	Load sterilizer per sterilization equipment operator instructions.
ii.	Select cycle: Standard
iii.	Operate machine per sterilization equipment operator instructions.

Functional testing was done by the manufacturer on units after 10 clean and sterilization cycles per the parameters in Table 1.

The sterilization parameters listed in Table 1 have been validated per ISO 14937:2009


Sterilization efficacy is the responsibility of the end user.

DISPOSAL

Recycle components or dispose of the product and its residual elements or waste items in accordance with the local laws and regulations.

DIRECTIONS FOR USE

- Use Map-iT® Diagnostic Mapping Catheter Extension Cable to connect an electrocardiogram recording system, pacing equipment or cardiac programmable stimulator and mapping catheter.
- Inspect the Map-iT® Diagnostic Mapping Catheter Extension Cable for damage (i.e. cuts, kinks, nicks, or crushed or elongated sections) prior to use. Do not use if damaged.
- Test the Map-iT® Diagnostic Mapping Catheter Extension Cable for leakage (lead to lead) and resistance (end-to-end) prior to use (leakage should not exceed 1.0 µA at 250 V and resistance should not exceed 10 Ω).
- To record intracardiac electrograms, connect the extension cable(s) to the catheter and then to the electrocardiogram recording system. Maintain the polarity of the proximal end connector pins of the Map-iT® Diagnostic Mapping Catheter Extension Cable(s) when connecting to the electrocardiogram recording system to comply with applicable safety standards.

- Consult the electrocardiogram recording system instructions for use for the proper connection of the patient grounding pad.
- Use care to isolate any unused connector pins. This will reduce the chances of developing accidental current pathways to the heart.

NOTES

NORMAL CONDITIONS OF USE

Operating temperature: -29 °C ~55 °C, Humidity rate: ≤ 85%.

This cable meets the requirements of IEC 60601-1 Ed. 3.1 for Type CF Applied Part.

Chamber dimensions, temperature distribution profile and the configuration of the products placed in the chamber are variables that must be considered in verifying the effectiveness of a specific sterilization cycle.

DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OR MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN, UNDER NO CIRCUMSTANCES SHALL ACCESS POINT TECHNOLOGIES EP INC., OR IT'S AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW.

Descriptions and specifications appearing in ACCESS POINT TECHNOLOGIES EP INC., printed matter, including this publication, are informational only and meant solely to generally describe the product at the time of manufacture and are not made or given as a warranty of the prescribed product in any way.

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Access Point Technologies EP Inc.