Map-iT® Steerable Ablation Catheter

INSTRUCTIONS FOR USE

Issue Date: 01/2026

I. Product Name: Map-iT® Steerable Ablation Catheter

II. Structural Performance: The radiofrequency ablation catheter is composed of a polyurethane catheter, platinum-iridium electrode tip, ring electrodes, push-pull handle, conductor wires, temperature sensor (thermocouple or thermistor), wire braiding, support tube, saline infusion tube (for the cold saline infusion RF ablation electrode catheter) and connector. Among those, the wires, temperature sensor (thermocouple or thermistor), wire braiding, saline infusion tube (for the cold saline infusion RF ablation electrode catheter) and support tube are buried in the catheter, while the temperature sensor is located in the center of the tip electrode to monitor the temperature of the tip electrode during the RF ablation. The catheter consists of a strong shaft and deflectable tube, which controls the bend of the catheter tip by tensing or relaxing the pull wires with a push rod connected to the push-pull handle to allow the catheter tip reach any part of the heart you wish to create a lesion. The catheter is connected to the RF generator via a connection cable with a corresponding interface, and RF energy is delivered where the catheter tip electrode comes in contact with the heart, to achieve ablation. There is a saline injection port with a standard Luer-lock interface connected to the proximal end of the catheter and it is used to inject saline through the lumen and out the distal end of the open end of saline infusion catheter, in order to flush the electrode to cool it. There are 6 infusion holes at the distal end of the tip electrode, uniformly distributed radially around the electrode tip. The end of the saline infusion lumen is connected to the infusion holes on the tip electrode. During the ablation application, the saline flows through the saline infusion lumen, and flows out from the electrode tip, and rinses and flushes the distal ablation electrode.

The catheter, for an **irrigated (CF) type** application, can be connected to a standard recording system and RF generator via connection cables with appropriate interfaces, and can operate only for a short time.

III. Specifications and Dimensions: For the specifications and dimensions of the Map-iT® Steerable Ablation Catheter, see Schedule 1.

IV. Product Usage: For radiofrequency ablation of heart tissue associated with tachycardias, or for temporary electrocardiography, mapping, pacing and stimulation.

V. Applicable Scope: The Map-iT® Steerable Ablation Catheter is applicable for cardiac electrophysiological mapping, stimulation and recording, and when used in conjunction with an RF generator, it can be used for intracardiac ablation, and treating tachycardias.


VI. Instructions:

The Map-iT® Steerable Ablation Catheter is a flexible electrode catheter that consists of a radiopaque polyurethane (PEBAX) insulator / wire braid, catheter equipped with platinum-iridium ring and tip electrodes, and electrode tip provided with a thermocouple or thermistor. A push-pull handle is installed at the proximal end of the catheter, and can be used to control the deflection of the distal end of the catheter.

The Map-iT® Steerable Ablation Catheter is used with an RF ablation generator, and the generator's technical specifications are as follows:

- Frequency: 450-575 kHz
- Voltage: 0-220 V
- Current: 0-2.0 A
- Power: 100-150 W (max.)
- Impedance cut-off limit: upper limit, 300Ω; lower limit: 25Ω.


The above output range refers to the equipment performance rather than the actual set values. The actual output set value is to be determined by the doctor based on experience and clinical requirements.

The temperature monitoring accuracy of the Map-iT® Steerable Ablation Catheter is $\pm 4.0\text{ }^{\circ}\text{C}$.

Connection between the Map-iT® Steerable Ablation Catheter and Radiofrequency Generator

1. Common Map-iT® Steerable Ablation Catheter

2. Saline infusion Map-iT® Steerable Ablation Catheter

The Map-iT® Steerable Ablation Catheter can be connected to various RF generators via connection cables. The catheter and connection cables are provided with quick plug connectors, and the part between the catheter and connection cable, as well as the part between the connection cable and RF generator are available for direct insertion and extraction. The following table shows the specifications of the catheter and connection cables selected according to the RF generator.

Model #	For use with what Map-IT® Ablation Catheter	RF Generator	Description	Length	Color
902402	Map-IT® Thermistor Ablation Catheter	EPT-1000	10 – 9-pin Extension Cable	300cm	Black
902404	Map-IT® Ablation Catheter	Cordis	10 – 10-pin Extension Cable	300cm	Black
902406	Map-IT® Ablation Catheter	IBI	10 – 14-pin Extension Cable	300cm	Black
902408	Map-IT® Ablation Catheter	Generic	10 – 7-pin Extension Cable	210cm	Black
902412	Map-IT® Thermocouple Ablation Catheter	Medtronic	10 – 10-pin Extension Cable	300cm	Black

Basic facilities and equipment required to be provided for the radiofrequency ablation procedure

- X-ray machine (DSA)
- Programmed electrical stimulator
- Intracardiac electrogram recording system
- Defibrillator
- RF generator
- RF ablation catheter
- Electrophysiological mapping catheter
- Connection cable
- Catheter introducer sheath
- Infusion pump
- 1ml or 2ml syringe

Before using this product, please read the following content carefully:

- Check the catheter's sterile packaging for intactness. Check the catheter's type of temperature sensor for its compatibility with the RF generator. Remove the catheter from the sterile package.
- If it is a saline infusion catheter, connect the saline infusion interface of the catheter to the infusion pump through an extension tube, and set the infusion pump speed to 60ml/min, start the infusion pump to flush the catheter infusion chamber, and observe the six infusion holes on the catheter electrode tip for a smooth flow and uniform infusion. If it does not do so, replace the catheter. When the air in the infusion chamber is drained, adjust the infusion pump speed to 2ml/min.
- A catheter sheath is used to insert the RF ablation catheter into the vasculature. Since the doctor needs to have easy control of the movement of the catheter within the heart chamber, it is recommended to use an 8F catheter sheath with a hemostasis valve for the use of a 7F radiofrequency ablation catheter.
- Upon insertion of the catheter, fluoroscopy shall always be used for guidance.
- Only the connecting cables for this electrode catheter can be used, and the RF ablation catheter is to be connected to the RF device or the electrode interface with a corresponding junction box.
- It is required to record the electrocardiogram, and the connection cable is to be connected to the electrode catheter. When the connection cable is connected to the ECG amplifier that meets the international safety standards, observe the polarity of the pins at the cable's proximal connector.
- If the device is to be used for temporary pacing, please connect the connecting cable to the electrode catheter. When the connection cable is connected to the temporary pacemaker in vitro that meets the international safety standards, observe the polarity of the pins at the cable's proximal connector.
- If this device is to be used for radiofrequency ablation without temperature monitoring, please connect the connection cable to the electrode catheter. When the connection cable is connected to the ECG amplifier that meets the international safety standards, observe the polarity of the pins at the cable's proximal connector. The doctor shall determine the exact location of the ablation point based on clinical electrophysiology experience, and help from X-ray guidance and mapping technologies. According to the above introductions in this manual, connect the proximal connector of the connection cable to the radiofrequency ablation instrument.
- When an RF ablation catheter with temperature monitoring is used, it is a must to select a radiofrequency ablation catheter based on the temperature monitoring characteristics of the RF generator used, and use the appropriate cables to connect the catheter to the radiofrequency ablation generator.
- For information about the proper use of the patient cable ground patches, please refer to the instructions for use of the radiofrequency ablation generator.

11. When the catheter is removed from the patient's body, it is required to keep the electrode catheter tip straight. Be careful to isolate all unused connector pins, which will reduce the chance of accidental current being introduced to the heart.
12. For manipulation of the electrode catheter tip, please push and pull the catheter handle knob. Forward pushing of the handle knob makes the catheter tip bend, while backward pulling of the handle knob makes the tip straighten.
13. The manipulation of the electrode catheter tip should always be performed under guidance with the fluoroscopy.
14. When the ablation site is determined, and the distal electrode has stable contact with the ablation site, the connection of the catheter tip electrode at the end shall be switched from the recording device to the RF generator, and the parameters of the RF device shall be set, and then the RF generator shall be started for radiofrequency ablation. If it is a saline infusion catheter, the infusion pump speed for ablation shall be adjusted to 17ml/min, and after an interruption of the RF current, the infusion pump shall be adjusted back to 2ml/min.
15. A catheter can be used to apply the RF current to the same or other sites once more. However, if the RF generator is turned off (due to the impedance or temperature), it is required to remove the catheter and clean any adhering clots on the electrode at the distal end, and then apply a second RF current. Sterile gauze soaked in sterile saline can be used to gently wipe the tip, and no forceful scrubbing or twisting of the electrode at the end is allowed, or it will damage the bonding of the electrode at the end and lead to its loosening. If it is a saline infusion catheter, it is required to confirm that the flush holes are not clogged before a second insertion. If the flushing hole is blocked:
 - a. Remove the catheter from the patient's body.
 - b. Inject some sterile saline solution into a 1ml or 2ml syringe, and connect it to the sidearm of the sheath.
 - c. Carefully inject the saline from the syringe into the catheter.
 - d. If necessary, repeat steps b) and c).
 - e. Once the flush holes have been completely cleaned, the catheter can be re-inserted into the patient's body. If the catheter is still blocked or does not work, do not use this catheter.
16. When the ablation is completed, pull the handle knob back to the bottom to straighten the tip, and then remove the catheter from the body.

VII. Warnings

-  The product should be stored in a cool and dry place.
-  Do not modify the device.
-  Please check the components before use.
-  It is not recommended for patients requiring long-term pacing.
-  Single-pole resistance should not exceed 50 ohms.
-  Do not use this device a second time. As it is impossible to completely remove the biological material and foreign matter, so the secondary use of the device may result in mechanical damage and / or a negative reaction in patients.
-  As catheter ablation can reach a certain intensity and persist for some time due to the X-ray beam if the fluorescence imaging, there is a potential danger to expose the patient and laboratory personnel to a large amount of X-rays, which may cause them to suffer acute radiation injury. Before carrying out the catheter ablation, pay full attention to the potential of radiation exposure related to the procedure, and take appropriate measures to minimize the exposure. Therefore, the use of this catheter in pregnant women must be carefully considered.
-  Patients undergoing atrioventricular node modulation or atrioventricular reentrant tachycardia ablation are vulnerable to dangers caused by inadvertent atrioventricular block. It is recommended to closely monitor the atrioventricular conduction during the RF power delivery. If partial or complete atrioventricular block is observed, immediately cut off the RF energy power.
-  A pacemaker or ICD may be affected by the negative impact of the RF signal. Be sure to remember:
 - a. Prepare a temporary pacing device in vitro for use during the ablation period;
 - b. Temporarily set the pacing system to the minimum output or the 000 mode, to minimize the risk of incorrect pacing;
 - c. While performing ablation, be particularly careful when it is very close to an atrial or ventricular permanent pacing lead;
 - d. After the ablation, a comprehensive pacing system analysis on all patients should be carried out;
 - e. While the power supply of the RF generator is on, the ICD should not be active.
-  During a retrograde aortic cannulation approach, fluoroscopy is required to ensure adequate visibility to avoid insertion of ablation catheter into the coronary arteries. Radiofrequency ablation with the catheter placed in a coronary artery will lead to a myocardial infarction and possible death.
-  Patients who consent to left atrial radiofrequency ablation should be closely monitored after the procedure for any clinical observation of an infarct phenomenon.

-  For patients who consent to cardiac ablation in the left atrium or atrial septum, it is recommended to administer anticoagulant therapy for the procedure; for patients who consent to right atrial procedures, the selective use of anticoagulant therapy during the procedure should be considered.
-  It is recommended that the RF power used in the catheter ablation should not exceed 30W.
-  The long-term risk of the extended exposure to fluoroscopy during prolonged procedures has not been determined. Therefore, the use of this catheter in children younger than adolescents must be carefully considered.
-  This device should only be used by physicians with adequate experience in intracardiac venous procedures, temporary pacing, and radiofrequency ablation techniques.
-  It is important carefully handle the catheter to avoid any cardiac damage, perforations, or tamponade. Insertion of the catheter must be conducted under fluoroscopic guidance. Do not push or pull the catheter if resistance is felt.
-  Do not insert or extract the catheter when the pullwire is tensioned, and the catheter tip is not straightened (confirm with fluoroscopy).
-  Observe the polarity.
-  Do not use it for direct current (DC) ablation.
-  It is important to prevent the connector from coming in touch with water or moisture, and only by doing so can you ensure the radiofrequency ablation electrode catheter functions normally.
-  The RF ablation electrode catheter should only be used with the appropriate connection cables, or it may result in injury to the patient or operator or equipment damage.
-  All RF ablation electrode catheters are equipped with a temperature sensor at the distal end. Appropriate connection cables are required. The end of the connection cable marked with "catheter" shall be connected to the RF ablation electrode catheter, and the end of the connection cable marked with "XX RF generator" should be connected to the "XX" RF generator.
-  The RF ablation electrode catheter is used with an RF generator. For detailed information about the operating instructions for the RF device, see the RF device operation manual.
-  As the different RF devices have their own specific temperature control mode, it is required to select a catheter based on the temperature control mode of the RF generator. The thermocouple-type RF device can only be equipped with a catheter marked with "TC / thermocouple"; thermistor RF device can only be equipped with a catheter marked with "THR / thermistor".
-  With irrigated ablation catheters the temperature sensor inside the electrode cannot reflect the electrode/tissue interface or tissue temperature due to the saline's cooling effect due to the flushing of the electrode. The temperature displayed on the RF generator refers to the cooled electrode temperature, rather than the tissue temperature. The temperature sensor is used to detect proper irrigation flow. Before applying RF current, a decrease in the electrode temperature confirms that the saline has already begun to flow through the ablation electrode. During the RF current application, recording the electrode temperature ensures the maintenance of the irrigation flow.
-  If performing left atrial ablation via a retrograde aortic approach, do not use a catheter 85cm in length.
-  Any sudden increase in the resistance during ablation usually indicates a conductor wire fracture and/ or distal coagulum formation. The catheter should be removed from the patient's body, and the distal section should be cleaned before continuing.
-  Use equipment isolated from the patients.
-  Do not insert the proximal connector pins of the patient cable or the electrode catheter directly into the electrical outlet (i.e. power sockets, wires, etc.).
-  You must use and maintain the RF ablation electrode catheter properly, and only this can ensure the device's normal electrical operation. When wiping off the catheter, pulling and/or twisting the device components may cause damage.
-  Never touch the distal electrode of the RF ablation electrode catheter when RF energy is being delivered from the RF generator, or it may harm the operator.
-  In essence, these are the general instructions. The doctor may have information that differs from the above summarized details basis on his/her personal clinical experience.
-  Be fully familiar with and master the operation of the push-pull control of the handle at the distal end.
-  Check for air bubbles before using saline irrigation during the procedure. Air bubbles in the saline irrigation may cause an embolism.

VIII. Potential complications:

Potential complications include, thromboembolisms, air embolisms, thrombophlebitis, catheter knotting, fractured wires, catheter's entangled in the chordae tendineae, pericardial tamponade, delayed or complete atrioventricular conduction block (AVB).

Preventive measures:

1. Flush the catheter with saline before putting it into the vasculature; high-risk patients should be treated prophylactically with heparin, and postoperatively with nadroparin calcium or warfarin.
2. Ensure a smooth venous backflow after the procedure.
3. Perform a pre-use check of the catheter integrity and stability, and that it can be removed from the vascular bifurcation or ascending aorta in preparation for the chance it could get caught.
4. Regulate the operation of the catheter.
5. When performing a transseptal puncture, do so under fluoroscopic guidance with a contrast medium, and use a transseptal puncture needle, dilator and outer sheath in the correct and a coordinated manner. Perform ultrasound guidance when necessary.
6. If III degree (complete) AVB occurs during the procedure, stop the ablation, and carry out an electrophysiological examination after one week.

IX. Contraindications:

1. The RF ablation electrode catheter is prohibited from use in compromised patients with a general infection in an active stage.
2. The atrial transseptal puncture procedure is prohibited from use in patients with a left atrial thrombus or myxoma and interatrial closure device or patch.
3. The aortic retrograde approach is prohibited from use in patients with a history of an aortic valve replacement surgery.
4. In patients that do not give consent to the cardiac catheterization procedure, the use of this product is prohibited.
5. If the patient suffers from an intracardiac mural thrombus or has undergone a ventricular dissection or atrial dissection within four weeks, this catheter is not allowed for use in radiofrequency ablation.

X. Transportation and storage:

1. Ambient temperature range: 0°C ~ 40°C;
2. Relative humidity range: 10% ~ 100%;
3. Atmospheric pressure range: 50kPa ~ 106 kPa.

XI. Service life: The Map-iT® Steerable Ablation Catheter uses ethylene oxide for sterilization, and its service life is for three years starting from the date of sterilization.

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Access Point Technologies EP Inc.