
# QuikCAS™ *with* Hemostasis Adapter IV:

Cardiodrive® *Sterile Components*

## INSTRUCTIONS FOR USE

**REF** 001-001100-1

|            |    |              |    |
|------------|----|--------------|----|
| English    | 5  | Italiano     | 40 |
| 中文         | 8  | Lietuvių k.  | 43 |
| 繁體中文       | 11 | Magyar       | 46 |
| 日本         | 14 | Nederlands   | 49 |
| 한국어        | 18 | Norsk        | 52 |
| Česky      | 21 | Polski       | 55 |
| Dansk      | 24 | Português    | 59 |
| Deutsch    | 27 | Русский язык | 62 |
| Español    | 31 | Română       | 66 |
| Eesti keel | 34 | Suomi        | 69 |
| Français   | 37 | Svenska      | 72 |

**R<sub>x</sub>**  
ONLY

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

DSP-0294 Rev: F

Effective Date: 08 28 2025

**SYMBOL LEGEND / 符号图例说明 / 符號圖例 / 凡例 / 기호 범례 / LEGENDA K SYMBOLŮM / SYMBOLFORKLARING / SYMBOLLEGENDE / LEYENDA DE SÍMBOLOS / SŪMBOLITE SELGITUS / LÉGENDE DES SYMBOLES / LEGENDA / SIMBOLIŲ APRAŠAS / JELMAGYARÁZAT / LEGENDA VAN SYMBOLEN / SYMBOLFORKLARING / OPIS SYMBOLI / LEGENDA DOS SÍMBOLOS / УСЛОВНЫЕ ОБОЗНАЧЕНИЯ / LEGENDA SIMBOLURILOR / MERKKIEN SELITYKSET / SYMBOLFÖRKLARING**

|                                                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                                                                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | <p>Catalog Number / 目录编号 / 目錄編號 / カタログ番号 / 카탈로그 번호 / Katalogové číslo / Katalognummer / Katalog-Nr. / Número de catálogo / Katalooginumber / Numéro de catalogue / Codice prodotto / Numeris kataloge / Katalógusszám / Catalogusnummer / Katalognummer / Numer katalogowy / Referência / Номер по каталогу / Număr de catalog / Luettelonumero / Katalognummer</p>                                                                                                                                                                                                                                                                     |    | <p>Sterilized using Ethylene Oxide / 已用环氧乙烷灭菌 / 使用環氧乙烷滅菌 / エチレンオキシド滅菌 / 에틸렌옥사이드로 멸균 처리됨 / Sterilizováno etylenoxidem / Steriliseret ved hjælp af ethylenoxid / Sterilisiert mit Ethylenoxid / Esterilizado con óxido de etileno / Steriliseeritud etüleenoksiidiga / Stérilisé à l'oxyde d'éthylène / Sterilizzato mediante ossido di etilene / Sterilizuota naudojant etileno oksidą / Etilén-oxiddal sterilizálva / Gesteriliseerd met behulp van ethyleenoxide / Steriliseret med etylenoksid / Wysterylizowano za pomocą tlenku etylenu / Esterilizado por óxido de etileno / Стерилизовано этиленоксидом / Sterilizat cu oxid de etilenă / Steriloitu etyleenioksidilla / Steriliserad med etylenoxid</p>     |
|    | <p>Lot Number / 产品批号 / 批次編號 / ロット番号 / 로트 번호 / Číslo šarže / Partinummer / Chargen-Nr. / Número de lote / Partii number / Numéro de lot / Numero di lotto / Partijos numeris / Gyártásítétel-szám / Partijnummer / Partinummer / Numer partii / Número de lote / См. руководство по эксплуатации / Номер партии / Număr de lot / Eränumero / Partinummer</p>                                                                                                                                                                                                                                                                             |    | <p>Double Sterile Barrier System / 双无菌屏障系统 / 雙層無菌隔絕系統 / ダブル無菌バリアシステム / 이중 무균 장벽 시스템 / Systém s dvojitou sterilní bariérou / Dobbelt steril barriresystem / System mit Doppel-Sterilbarriere / Sistema de barrera estéril doble / Kahekordse steriilse kaitse süsteem / Système à double barrière stérile / Sistema a doppia barriera sterile / Dviguba sterilaus barjero 2istema / Dupla steril gátrendszer / Systeem met dubbele steriele barrière / Dobbelt steril barriresystem / System podwójnej bariery sterylnej / Sistema de barreira estéril dupla / Система двойного стерильного барьера / Sistem de barieră dublu sterilă / Kaksinkertainen steriili estojärjestelmä / Dubbelt steril barriärsystem</p> |
|  | <p>Consult Operating Instructions / 请查阅操作说明 / 參閱操作說明 / 取扱説明書を参照 / 작동 지침 참조 / Nahlédněte do návodu k použití / Se brugsanvisning / Bedienungsanleitung konsultieren / Consultar las instrucciones de funcionamiento / Lugege kasutusjuhendit / Consulter le mode d'emploi / Consultare le istruzioni per il funzionamento / Žr. naudojimo instrukciją / Ellenőrizze az üzemeltetési utasításokat / Raadpleeg de bedieningsinstructies / Les bruksanvisningen / Zapoznać się z instrukcją obsługi / Consultar as Instruções de Funcionamento / Consultați instrucțiunile de utilizare / Tutustu käyttöohjeisiin / Se bruksanvisningen</p> |  | <p>Single Use Only / 仅限一次性使用 / 僅供單次使用 / 使い捨てのみ / 일회용 / Pouze na jedno použití / Må ikke resteriliseres / Kun til engangsbrug / Verwendbar bis / Nur zum einmaligen Gebrauch / Para un solo uso / Ainult ühekordseks kasutamiseks / Réserve à un usage unique / Esclusivamente monouso / Tik vienkartiniam naudojimui / Kizárolag egyszer használatos / Uitsluitend voor eenmalig gebruik / Kun til engangsbruk / Do jednorazowego użytku / Apenas para uma única utilização / Только для одноразового применения / De unică folosință / Ei saa käyttää uudelleen / Endast för engångsbruk</p>                                                                                                                      |
|  | <p>Use By / 使用截止日期 / 使用期限 / 使用者 / 유효 기한 / Použitelné do / Anvendes inden / Verwendbar bis / Fecha de caducidad / Kasutada enne / Utiliser avant le / Data di scadenza / Sunaudoti iki / Szavatossági idő / Uiterste gebruiksdatum / Brukes innen / Termin ważności / Utilizar até / Исползовать до / Expiră la / Viimeinen käyttöpäivä / Använd före</p>                                                                                                                                                                                                                                                                                |  | <p>Do Not Resterilize / 请勿重新消毒 / 請勿重複滅菌 / 再滅菌しないこと / 재멸균하지 마십시오 / Opakovaně nesterilizujte / Må ikke resteriliseres / Nicht erneut sterilisieren / No 2teril a esterilizar / Árge steriliseerige korduvalt / Ne pas restériliser / Non risterilizzare / Nesterilizuoti pakartotina / Újrasterilizálása tilos / Niet opnieuw steriliseren / Ikke 2terilizer på nytt / Nie sterylizować ponownie / Não reesterilizar / He sterilizovat повторно / A nu se reesteriliza / Ei saa steriloida uudelleen / Får ej omsteriliseras</p>                                                                                                                                                                                     |

|                                                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |                                                                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|     | <p>EU Conformity Assessed / 经过欧盟符合性评估 / 歐盟合格認證 / 英国適合性評價済み / EU 적합성 평가 완료 / Posouzena shoda v EU / EU-overensstemmelsesvurderet / EU-Konformität geprüft / Evaluación de conformidad de la UE / ELI vastavushindamine / Conformité de l'UE évaluée / Conformità UE valutata / Įvertinta ES atitiktis / EU megfelelőségi tanúsítvány / EU-conformiteitsbeoordeling / EU-samsvar vurdert / Potwierdzona zgodność z przepisami obowiązującymi w Unii Europejskiej / Conformidade da UE avaliada / Оценка соответствия требованиям ЕС пройдена / Evaluat pentru conformitate în UE / EU:n vaatimustenmukaisuus arvioitu / EU-överensstämmelse bedömd</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |    | <p>Packaging Unit / 包装单位 / 包裝單元 / 包装ユニット / 포장 단위 / Balení / Emballeringsenhed / Verpackungseinheit / Unidad de embalaje / Pakkeühik / Unité d'emballage / Unità di confezionamento / Pakuotė / Csomagolási egység / Verpakkingsseenheid / Emballasjeenhet / Jednostka opakowania / Unidade de embalagem / Упаковочная единица / Unitate de ambalare / Pakkausyksikkö / Förpackning</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|    | <p>Do Not Use if Package is Damaged or Opened / 如果包装打开或损坏 · 请勿使用 / 若包裝已打開或有毀損，請勿使用 / パッケージが破損または開封されている場合は使用しない / 포장이 손상되었거나 개봉된 경우 사용하지 마십시오 / Pokud je obal otevřený nebo poškozený, zařízení nepoužívejte / Må ikke anvendes, hvis pakken er beskadiget eller åbnet / Bei beschädigter oder geöffneter Verpackung nicht verwenden / No 3tilizer si el embalaje está dañado o abierto / Mitte kasutada, kui pakend on kahjustatud või avatud / Ne pas utiliser si l'emballage est endommagé ou ouvert / Non utilizzare se la confezione è danneggiata o aperta / Nenaudoti, jeigu sugadinta arba atidaryta pakuotė / Tilos felhasználni, ha a csomagolás fel van bontva vagy sérült / Niet gebruiken als verpakking beschadigd of geopend is / Må ikke brukes hvis pakningen er skadet eller åpnet / Nie używać, jeżeli opakowanie jest uszkodzone lub otwarte / Não utilizar se a embalagem estiver danificada ou aberta / Не использовать, если упаковка повреждена или вскрыта / A nu se utiliza dacă ambalajul este deteriorat sau a fost deschis / Ei saa käyttää, jos pakkaus on vaurioitunut tai avattu / Får ej användas om förpackningen är skadad eller öppnad</p> |    | <p>Importer / 进口商 / 進口商 / 輸入業者 / 수입업체 / Dovožce / Importør / Importeur / Importador / Importija / Importateur / Importatore / Importuotojas / Importör / Importeur / Importør / Importer / Importador / Импортёр / Evaluat pentru conformitate în Regatul Unit / Importator / Maahantuoja / Importör</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|   | <p>Activation Symbol / 激活符号 / 啟動符號 / アクチベーションシンボル / 활성화 기호 / Aktivací symbol / Aktiveringssymbol / Aktivierungssymbol / Símbolo de activación / Aktiveerimissymbol / Symbole d'activation / Simbolo di attivazione / Aktyvinimo simbolis / Aktiválási szimbólum / Activeringssymbol / Aktiveringssymbol / Symbol aktywacji / Símbolo de ativação / Символ активации / Simbol de activare / Aktivointisymboli / Aktiveringssymbol</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |  | <p>Medical Device / 医疗设备 / 醫療裝置 / 医療機器 / 의료 기기 / Zdravotnický prostředek / Medicinsk udstyr / Medizinisches Gerät / Dispositivo médico / Meditsiiniseade / Dispositif médical / Dispositivo medico / Medicinos prietaisas / Orvostechnikai eszköz / Medisch apparaat / Medisinsk utstyr / Wyrób medyczny / Dispositivo médico / Медицинское изделие / Dispozitiv medical / Lääkinnälinen laite / Medicinsk enhet</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|  | <p>Keep Dry / 保持干燥 / 保持乾燥 / 水濡厳禁 / 건조한 장소에 보관하십시오 / Uchovávejte v suchu / Opbevares tørt / Trocken lagern / Manter seco / Hoida kuivas / Conserver au sec / Manterene asciutto / Laikyti sausoje vietoje / Szárazon tartandó / Droog houden / Må holdes tørr / Utrzymywać produkt w stanie suchym / Manter seco / Хранить в сухом месте / A se păstra la loc uscat / Säilytettävä kuivana / Förvaras torrt</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |  | <p>Authorized Representative in the European Community / 欧盟授权代表 / 歐盟授權代表 / 欧州委任代理人 / 유럽공동체(EC) 공식 대리점 / Autorizovaný zástupce v Evropském společenství / Autoriseret repræsentant i EU / Autorisierte Vertretung in der Europäischen Gemeinschaft / Representante autorizado en la Comunidad Europea / Volitatus esindaja Euroopa Ühenduses / Représentant agréé établi dans la Communauté européenne / Rappresentante autorizzato per la Comunità Europea / Įgaliotasis atstovas Europos Bendrijoje / Meghatalmazott képviselő az Európai Közösségben / Gemachtigd vertegenwoordiger in de Europese Gemeenschap / Autoriseret representant i EU / Autoryzowany przedstawiciel we Wspólnocie Europejskiej / Representante autorizado na Comunidade Europeia / Уполномоченный представитель в Европейском сообществе / Reprezentant autorizat în Comunitatea Europeană / Valtuutettu</p> |

|                                                                                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                                                                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                   |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |                                                                                    | edustaja Euroopan yhteisössä / Auktoriserad EU-representant                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|  | <p>Keep Away from Sunlight / 避免日照 / / 避免日光照射 / 日光を避ける / 직사광선에 노출하지 마십시오 / Nevystavujte slunečnímu světlu / Opbevares beskyttet mod sollys / Vor Sonnenlicht schützen / Mantener alejado de la luz solar / Hoida eemal päikesevalgusest / Conserver à l'abri de la lumière du soleil / Tenere lontano dalla luce solare / Laikyti atokiai nuo tiesioginės saulės šviesos / Napfénytől védve tartandó / Niet blootstellen aan zonlicht / Må beskyttes mot sollys / Przechowywać z dala od promieni słonecznych / Manter afastado da luz solar / Избегайте попадания прямых солнечных лучей / A se feri de lumina soarelui / Suojattava auringonvalolta / Skyddas från solljus</p>                                                                                                                                                                                                                                                                                                                 |  | <p>Prescription Device Only / 仅限处方设备 / 本設備僅限遵醫囑使用 / 処方医療機器のみ / 처방 장치 전용 / Prostředek pouze na lékařský předpis / Kun receptpligtigt udstyr / Verschreibungspflichtiges Gerät / Dispositivo de venta bajo receta / Ainult retseptiga seade / Dispositif sur ordonnance uniquement / Dispositivo disponibile solo su prescrizione / Prietaisas naudojamas tik pagal paskyrimą / Kizárolag orvosi rendeltvényre / Hulpmiddel is uitsluitend verkrijgbaar op recept / Reseptbelagt enhet / Wyrób tylko na receptę / Dispositivo apenas com receita médica / Изделие может быть приобретено только по назначению врача / Dispozitiv eliberat numai pe bază de prescripție medicală / Vain reseptillä määrättävä laite / Produkt endast mot recept</p> |
|  | <p>Protect from Heat and Radioactive Sources / 远离热源和放射源 / 遠離熱源及放射源 / 熱および放射線源から保護する / 열원 및 방사능원으로부터 보호하십시오 / Chraňte před zdroji tepla a radioaktivního záření / Beskyttes mod varme og 4radioactive kilder / Vor Hitze und radioaktiven Strahlungsquellen schützen / Proteger del calor y de fuentes radiactivas / Kaitsta kuumuse ja radioaktiivse kiirguse allikate eest / Protéger de la chaleur et des sources radioactives / Proteggere dal calore e da fonti radioattive / Saugoti nuo karščio ir radioaktyviųjų šaltinių / Övja a hőtől és radioaktív sugárforrásoktól / Bescherm tegen warmtebronnen en radioactieve bronnen / Må beskyttes mot varme og 4radioactive kilder / Chronić przed źródłami ciepła i promieniowania / Proteger do calor e fontes radioativas / Храните вдали от источников тепла и радиоактивного излучения / A se feri de căldură și surse de radioactivitate / Suojattava lämmöltä ja radioaktiivisuudelta / Skyddas från värme och radioaktiva källor</p> |  | <p>Manufacturer / 制造商 / 製造商 / メーカー / 제조업체 / Výrobce / Producent / Hersteller / Fabricante / Tootja / Fabricant / Produttore / Gamintojas / Gyártó / Merk / Produsent / Producent / Fabricante / Производитель / Producător / Valmistaja / Tillverkare</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |

## ENGLISH

### DEVICE DESCRIPTION

*QuikCAS* is a part of the *Cardiodrive* automated catheter advancement system (CAS). *QuikCAS* advances or retracts a compatible magnetic electrophysiology (EP) catheter when used in conjunction with a Stereotaxis Robotic Magnetic Navigation System (RMNS).


The Hemostasis Adapter is attached to a hemostasis introducer, which an EP catheter is fed through. The Advancer Unit is attached to the *Cardiodrive* motor. The EP catheter is controlled remotely via a user interface (UI) located at the patient table, or in the control room.

### STERILE DISPOSABLE COMPONENTS:

Hemostasis Adapter

Pre-Assembled Advancer Unit

### ADDITIONAL EQUIPMENT

The *Cardiodrive* hardware components and Stereotaxis RMNS should be used with *QuikCAS*. For more information, refer to the applicable RMNS User Guide.

### HOW SUPPLIED

Packaging is designed to maintain sterility according to the expiration date on the label. Do not use if the expiration date has passed, the package is opened or damaged, or if the labeling is illegible.

### HANDLING AND STORAGE

Store at ambient temperatures in a dry and dark location.

### COMPATIBLE DEVICES

The following disposable items are specified for *QuikCAS* operation, and sold separately.

- Compatible magnetic EP catheters:
  - NAVISTAR® RMT Catheter (7F)
  - NAVISTAR® RMT THERMOCOOL® Catheter (8F)
  - CELSIUS® RMT Catheter (7F)
  - CELSIUS® RMT THERMOCOOL® Catheter (8F)
  - MedFact MagnoFlush Catheter (8F)
  - Magnetic Interventional Ablation Catheter (MAGiC™)
- Designated introducer sheaths or guiding sheaths:
  - Abbott Fast-Cath™ 8.5F Hemostasis Introducer (P/N: 406123, 406976)
  - Abbott Fast-Cath 8.5F Transseptal Guiding Introducer (P/N: 406804, 406902)
  - Abbott Swartz™ 8.5F Braided Transseptal Guiding Introducer (P/N: 407451, 407452, 407453, and 407454)

- Abbott Agilis™ NxT 8.5F Steerable Introducer (P/N: G402047, G408318, G408319, G408332)
- Biosense Webster MOBICATH® 8.5F Guiding Sheath Small Curve (P/N: D140010, D140011)
- Biosense Webster CARTO VIZIGO® 8.5F Bi-Directional Guidance Sheath (P/N: D138501, D138502, and D138503)
- Teleflex Super Arrow-Flex® 9F Sheath Introducer (P/N: CL-07965 and CL-07980)

### INDICATIONS

In the United States and Canada, the following Indications are applicable for the *Cardiodrive* system.

- The Stereotaxis *Cardiodrive* automated catheter advancement system (CAS) is intended for automatically advancing and retracting only compatible magnetic electrophysiology (EP) mapping and ablation catheters inside the patient's heart when used in conjunction with a Stereotaxis Robotic Magnetic Navigation System.
- The *Cardiodrive* system is not intended to advance the EP mapping and ablation catheters through the coronary vasculature or the coronary sinus.

In all other geographies, the following Indications are applicable for the *Cardiodrive* system.

- The Stereotaxis *Cardiodrive* automated catheter advancement system (CAS) is intended for automatically advancing and retracting only compatible magnetic electrophysiology (EP) mapping and ablation catheters inside the patient's heart and pericardial space when used in conjunction with a Stereotaxis Robotic Magnetic Navigation System.

### INTENDED PATIENT POPULATION

The intended patient population for the *QuikCAS* device is patients undergoing EP mapping and ablation procedures of the heart.

### INTENDED USERS

The *QuikCAS* device should be used only by qualified medical professionals who have been thoroughly trained in its use.

### WARNINGS

- The *QuikCAS* sterile components are intended for **SINGLE USE ONLY**. Do not reuse, reprocess, or resterilize. Reuse can compromise the performance characteristics and may result in patient infection.
- The components of the *QuikCAS* are sterile and should be handled under sterile conditions. Verify that the package integrity has been maintained so that the sterility of the device has not been compromised.
- The *QuikCAS* system should only be operated for advancement or retraction while the catheter is under direct fluoroscopic visualization.
- Ensure all blood is flushed out of the Hemostasis Adapter during the patient setup.
- Inadvertent advancement or retraction of the *QuikCAS* may cause contact with sensitive heart tissue leading to a possible arrhythmia.
- After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.

### PRECAUTIONS

- The *QuikCAS* system should be used only by properly trained physicians.
- Always reduce the magnetic fields before attempting to connect/disconnect the *QuikCAS* components or to insert/remove the catheter.
- The RMNS, used in conjunction with *QuikCAS*, produces a strong magnetic field during operation. This magnetic field may cause metal objects to become airborne projectiles if such objects are not kept sufficiently isolated from the magnetic field.
- Refer to the RMNS User Guide for user precautions.
- Restrain the patient's leg to avoid excess bending, which may interfere with *QuikCAS* performance and place excessive pressure on the vascular access site.
- Carefully examine the device for defects and verify proper device function and integrity before catheter insertion into the patient.

- Do not use any *QuikCAS* components that have been damaged.
- Use the sterile components of the *QuikCAS* device before the expiration date printed on the product package.

## PREPARATION FOR USE

1. Follow instructions in the *Cardiodrive* Installation, Operation, and Service Manual for setup of the User Interface and hardware components of the system.
2. Remove the *QuikCAS* Sterile Components tray from the box.
3. Navigant™ software version 3.0.6.X and later will require an activation bar code to be scanned. This bar code can be found on the box or tray label located by the “Activation Code” symbol (see below). When prompted by *Navigant*, use the bar code scanner to read the bar code.

4. Using sterile technique, peel back the lid from the outer tray, remove the inner tray from the outer tray, and peel back the lid from the inner tray.
5. Inspect the parts to confirm they are undamaged. If damaged, use a new *QuikCAS* for the procedure.
6. Remove the pre-assembled Advancer Unit from the tray. Orient the Advancer Unit so the Flexible Drive Cable (on the Advancer Unit) is pointing away from the intended introducer access site.
7. Position the Advancer Unit on the patient's leg to optimize horizontal placement.
8. Pull back the lever arm on the Advancer Unit to put it in the UNLOCKED position.


**CAUTION:** Sterile to non-sterile area connection required for next steps.

9. Attach the Flexible Drive Cable to the Motor Assembly.
  - a. Hand the end of the sterile cable to personnel **outside the sterile field**.
  - b. Personnel **outside the sterile field** must insert the cable through the Luer-Lock fitting located on the Motor Assembly. Ensure the cable aligns with the square profile of the coupler on the Motor Assembly shaft.
  - c. Tighten the Luer-Lock fitting to secure the connection.
10. Position the Motor Assembly near the foot of the bed such that the Flexible Drive Cable is in an “S” shape on the patient drape. Reduce the slack in cable as much as possible without pulling against the Advancer Unit.
11. Proceed to Patient Setup.

## PATIENT SETUP

**WARNING:** Immobilize or restrain the patient's leg to avoid excess bending and excessive pressure on the vascular access site.

1. The Hemostasis Adapter is designed for use with different commercially available introducer sheaths. Prepare and advance the introducer in accordance with the manufacturer's Instructions For Use.
  - a. The following images detail the introducer sheaths and guiding sheaths and the appropriate slots
2. Remove the Hemostasis Adapter from the sterile tray.
3. Prepare the compatible magnetic EP catheter in accordance with the manufacturer's Instructions For Use.
4. Insert the magnetic EP catheter past the tip of the Hemostasis Adapter.
5. Insert the non flexible tip of the magnetic EP catheter into the hemostasis valve of the sheath.
6. Position the Hemostasis Adapter so the compatible slot is available for easy introducer insertion.
7. Twist the introducer into the correct slot to lock the sidearm into the slot. If there is resistance, the components are not correctly aligned.

2. Abbott *Fast-Cath* 8.5F Hemostasis Introducer  
Abbott *Fast-Cath* 8.5F Transseptal Guiding Introducer  
Teleflex *Super Arrow-Flex* 9F Sheath Introducer


3. Biosense Webster *MOBICATH* 8.5F Guiding Sheath  
Biosense Webster *CARTO VIZIGO* 8.5F Bi-directional Guiding Sheath

8. Gently move the Hemostasis Adapter handles all the way forward. The front tube should extend into the compatible introducer's hemostasis valve. The introducer's hemostasis valve is now open.

**CAUTION:** Use the Hemostasis Adapter to **open** the hemostasis valve before continuing to pass the compatible magnetic EP catheter through the hemostasis valve.


9. Slowly advance the compatible magnetic EP catheter through the Hemostasis Adapter until the flexible section is completely inside the introducer. If there is resistance, ensure the handles are pushed all the way forward.
10. Grab the handles of the Hemostasis Adapter and move them backward/down until the cap on the end of the Hemostasis Adapter protrudes out of the adapter.

**NOTE:** In this position, the front tube is no longer extending into the introducer and the hemostasis valve is now closed around the catheter.

11. Attach a syringe to the sidearm of the introducer.
12. Open up the sidearm valve and **SLOWLY** drawback blood into the syringe until all the air is out of the system.
13. Close the sidearm valve and disconnect the syringe from the sidearm. Evacuate the syringe and fill the syringe with sterile saline. Remove all air from the saline-filled syringe and reattach the syringe to the sidearm. Open the sidearm, and flush the introducer.
14. Using the proximal end of the catheter, continue to feed it through the introducer until the catheter tip is slightly protruding beyond the end of the introducer.
15. The cap should be protruding out of the hemostasis adapter and ready to be seated into the Advancer unit. If the cap is not protruding, pull the Hemostasis Adapter handles all the way back.
16. Lower the catheter and the Hemostasis Adapter into the Advancer Unit:
  - a. Seat the proximal end of the Hemostasis Adapter into the slot of the Advancer Unit. The slot is at the end of the catheter channel, on the opposite side of the drive cable.
  - b. Position the catheter into the bottom of the catheter channel. If the catheter is not fully seated in the channel, it may be pinched when the Advancer Unit is closed.

1. Abbott *Swartz* 8.5F Braided Transseptal Guiding Introducer  
Abbott *Agilis* NxT 8.5F Steerable Introducer

17. Close the Advancer Unit. Verify the catheter is fully seated in the catheter channel and the lever is in the LOCKED position.
18. Loosen the thumbscrews on the Advancer Unit to allow vertical adjustment while positioning onto the patient.
19. Adjust Advancer Unit position to allow the Hemostasis Adapter tubing to be as straight as possible. Retighten the thumbscrews.
20. Ensure the Flexible Drive Cable is in an "S" shape on the patient drape. Secure the Flexible Drive Cable to the drape with tape.
21. If the catheter needs to be advanced further to reach the heart, move the lever arm to the UNLOCKED position on the Advancer Unit and manually advance the catheter until the tip has entered the heart.

**CAUTION:** Do not use the *QuikCAS* system to advance the catheter tip through the introducer. The force needed to advance the catheter tip may be more than the *QuikCAS* system is designed to deliver.

22. Move the lever arm to the LOCKED position on the Advancer Unit. Verify that the catheter is fully seated in the catheter channel and the lever is in the LOCKED position. This connects the catheter to the Advancer Unit and enables the use of the *QuikCAS*. Ensure the catheter is not pinched in the Advancer Unit.
23. Proceed to *QuikCAS* Operation.

**CAUTION:** Never attempt to move the Hemostasis handles forward while attached to the *QuikCAS*. This will damage the Hemostasis Adapter.

#### QUIKCAS OPERATION

1. With direction from the physician, the catheter may be advanced or retracted. Refer to the *Cardiodrive* User Guide for instruction of the entire system.
2. Aspirate and back flush the introducer during the EP procedure in accordance with the manufacturer's Instructions for Use.
3. Follow standard practice of using a continuous drip of anticoagulant fluid under pressure through the introducer's sidearm, in accordance with the manufacturer's Instructions for Use.

#### TROUBLESHOOTING

1. If the catheter ceases to advance or retract smoothly:
  - a. Verify the catheter tip has not been retracted into the introducer distal curved section. If so, open the Advancer Unit and manually advance the catheter until the tip has entered the heart.
  - b. Verify the Hemostasis Adapter and the catheter channel of the Advancer Unit are still aligned with the introducer. Loosen the thumbscrew and adjust as necessary.
  - c. Aspirate the introducer to remove fibrin buildup, in accordance with the manufacturer's Instructions for Use.
  - d. Observe the inner drive shaft of the Flexible Drive Cable when an advance or retract is commanded. If the inner drive shaft is moving, go to the next step. If the inner drive shaft is not moving, refer to the *Cardiodrive* User Guide.
  - e. Move the lever to the UNLOCKED position on the Advancer Unit. Slowly rotate the shaft of the compatible magnetic EP catheter approximately 90 degrees. Close the Advancer Unit. If the catheter has slipped in the Advancer Unit, a "smooth" spot may be created on the catheter shaft. Rotating the catheter shaft will be necessary to avoid this "smooth" spot during the balance of the EP procedure.
2. If the activation barcode does not activate the system:
  - a. Verify that the expiration date has not passed. If the device has expired, use a new device.
  - b. Verify that the device has not already been used. The device is intended for single use only.
  - c. Verify that the activation bar code has not been damaged. If it has been damaged, type the numerical code into *Navigant*.

3. If the catheter needs to be removed (during or after the procedure), remove the Hemostasis Adapter from the *QuikCAS*. Push the handles of the Hemostasis Adapter forward to insert the front tube into the hemostasis valve of the introducer. From the proximal end of the Adapter, gently pull the catheter out of the Adapter and Introducer.

**NOTE:** If *QuikCAS* fails to operate correctly after following these procedures, call Stereotaxis TST for further assistance.

#### STEREOTAXIS TECHNICAL SUPPORT

For technical support, please contact Stereotaxis TeleRobotic Support Team (TST) at 1-866-269-5268 or 1-314-678-6200 or email [tst@stereotaxis.com](mailto:tst@stereotaxis.com).

#### NOTICE TO THE USER AND/OR PATIENT

Any adverse event or malfunction that occurred in relation to the device should be reported to the manufacturer. Serious adverse events should be reported to the competent authority of the Member State in which the user and/or patient is established.

#### TRADEMARKS

Stereotaxis, the Stereotaxis logo, *QuikCAS*, *Cardiodrive*, and *Navigant* are registered trademarks of Stereotaxis, Inc. in the USA and other countries. All other brand names, product names, or trademarks are the property of their respective owners.

#### PATENTS

Manufactured under one or more of the following patents:  
 United States: 7,276,044; 7,635,342; 7,766,856  
 European Patents: EP 1 389 958 (B1); EP 1 781 364 (A2)

#### DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

STEREOTAXIS MAKES NO WARRANTIES, EXPRESS OR IMPLIED, WITH REGARD TO THE DEVICE DESCRIBED IN THIS DOCUMENT. STEREOTAXIS DISCLAIMS ALL REPRESENTATIONS OR WARRANTIES, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR USE, TITLE OR NON-INFRINGEMENT, ARISING BY STATUTE OR IN LAW, OR ARISING FROM A COURSE OF CONDUCT, COURSE OF DEALING OR USAGE OF TRADE.

STEREOTAXIS, INC. SHALL NOT BE LIABLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW.

STEROTAXIS IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES RESULTING FROM REUSE OF THIS DEVICE.

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© Stereotaxis, Inc. 2025
GmbH)