

HIS RV™ DIAGNOSTIC CATHETER

Multi-Electrode Catheter for Simultaneous HIS & RV Recording



Intentional Electrode Placement

Receive high quality electrogram data from four distal RV and six proximal HIS electrodes.



Streamlined Procedure

With a single catheter, pace the RV and record the AV node activation for a more efficient procedure.




Enhanced Access & Patient Safety

Reduce load on vasculature and preserve the opportunity for additional catheter access.


The Map-iT HIS RV, also known as the HIS-RV Diagnostic Catheter, allows physicians to use one catheter and a single access point for both HIS and RV recordings, simplifying procedures and potentially reducing the risk of HV interval prolongation. Map-iT catheters are designed with a braid-reinforced shaft for high pushability, easy placement, and consistent handling, even during long procedures.

Specifications

| | |
|-----------------------------|----------------------|
| French Size | 5F or 6F |
| Length | 120cm |
| Electrodes | 4 Distal, 6 Proximal |
| HIS to RV Electrode Spacing | 20, 23, 24, 29mm |
| Curve Type | MPA |



Superior Shaft Construction
High-performance torque control, pushability, and durability



4 Distal RV Electrodes

6 Proximal HIS Electrodes

Smooth shaft electrode interface and platinum iridium electrodes for high quality signals

5F Decapolar Fixed Diagnostic Catheter

| | |
|----------------|-------|
| French | 5F |
| Electrodes | 10 |
| Ring Electrode | 1mm |
| Tip Electrode | 2mm |
| Curve | MPA |
| Usable Length | 120cm |

| Reference No. | Electrode Spacing (mm) |
|---------------|------------------------|
| 901579 | 4-4-4-23-4-4-4-4-4 |
| 901580 | 4-4-4-29-4-4-4-4-4 |
| 901583 | 5-5-5-20-3-3-3-3-3 |
| 901584 | 5-5-5-24-3-3-3-3-3 |

Connection Cable Reference No. 901957

Made in the USA

6F Decapolar Fixed Diagnostic Catheter

| | |
|----------------|-------|
| French | 6F |
| Electrodes | 10 |
| Ring Electrode | 1mm |
| Tip Electrode | 2mm |
| Curve | MPA |
| Usable Length | 120cm |

| Reference No. | Electrode Spacing (mm) |
|---------------|------------------------|
| 901823 | 4-4-4-23-4-4-4-4-4 |
| 901824 | 4-4-4-29-4-4-4-4-4 |
| 901827 | 5-5-5-20-3-3-3-3-3 |
| 901828 | 5-5-5-24-3-3-3-3-3 |

Connection Cable Reference No. 901957

Made in the USA



To Order

✉ Map-iT@Stereotaxis.com ☎ +1 (763) 432-9004

Follow Us

✕ in EP Map-iT