



# CS REACH™

## DIAGNOSTIC CATHETER

Coronary Sinus Catheter with 3D Curve



### Access Challenging Anatomy

Designed for easy and rapid access through challenging coronary sinus anatomy: CS ostium displacement, Thebesian valve, Chiari network



### 90 Degree Curve

Multi-plane curve allows 90° approach to the CS ostium




### Stable Anchoring

3D curve allows for anchoring and maintains handling characteristics throughout long procedures

The Map-iT CS Reach, also known as the CS-Reach Diagnostic Catheter, is designed with a 90-degree curve for effortless navigation through demanding coronary sinus anatomy to the CS ostium. Once in place, this 3D curve provides stable anchoring for the duration of the procedure. All Map-iT catheters are designed with a braid-reinforced shaft for high pushability, easy placement, and consistent handling, even during long procedures.

## Specifications

|                   |                  |
|-------------------|------------------|
| French Size       | 5F or 6F         |
| Length            | 110cm            |
| Electrodes        | 10               |
| Electrode Spacing | 2-5-2mm, 2-8-2mm |
| Curve Type        | MPD-SL           |
| Deflection        | Multi-plane      |



## 5F & 6F Decapolar Steerable Diagnostic Catheter

|                |        |
|----------------|--------|
| Electrodes     | 10     |
| Ring Electrode | 1mm    |
| Tip Electrode  | 2mm    |
| Curve          | MPD-SL |
| Usable Length  | 110cm  |

| Reference No. | Electrode Spacing (mm) | French |
|---------------|------------------------|--------|
| 201707        | 2-5-2                  | 5      |
| 201708        | 2-8-2                  | 5      |
| 201717        | 2-5-2                  | 6      |
| 201718        | 2-8-2                  | 6      |

Connection Cable Reference No. 901957

Made in the USA



### To Order

✉ Map-iT@Stereotaxis.com ☎ +1 (763) 432-9004

### Follow Us

✕ in EP Map-iT